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
The VITEK® 2 Gram Negative Susceptibility Card is intended for use with the VITEK 2 Systems in clinical laboratories as an *in vitro* test to determine the susceptibility of clinically significant aerobic gram-negative bacilli to antimicrobial agents when used as instructed in the Product Information manual.

Applicable warranty is fully set forth in the instructions for use.

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.
**N/A = No specific FDA Indications for Use available
① ② etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in travelers diarrhea)

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AST-N270 REF 414 162

LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):

- Amikacin: *Acinetobacter baumannii*
- Amoxicillin/Clavulanic Acid: *Providencia* spp.
- Ampicillin: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*

Perform an alternative method of testing prior to reporting of results when an intermediate or resistant result is obtained with the following antibiotic/organism combination(s):

- Amoxicillin/Clavulanic Acid: *Proteus* spp.

The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:

- Ertapenem: *E. coli*, *K. pneumoniae*, *C. freundii*, *C. koseri*, *E. aerogenes*, *E. cloacae*, *K. oxytoca* (excluding ESBL-producing strains), *M. morgani*, *P. mirabilis*, *P. vulgaris*, *S. marcescens*

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NOTE: A result for an antibiotic/organism combination which may have a limitation may be suppressed from reporting. This can be accomplished through use of Conditional Antibiotic Reporting (CAR) rules. Instructions for use are contained in the bioLiaison® Procedures Manual under the section titled "Using Conditional Antibiotic Reporting".

↓ MUST enter Package Insert Barcode into “Flex Panel Entry” program before first use of this Susceptibility Card.

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Vitek and bioLiaison are used, pending and/or registered trademarks of bioMérieux. ATCC is a used, pending and/or registered trademark of American Type Culture Collection. CLSI is a used, pending and/or registered trademark of Clinical and Laboratory Standards Institute, Inc.

Product covered by one or more of U.S. Patent Nos. D414,272; D437,797; 5,609,828; 5,746,980 ; 5,804,437 ; 5,869,005; 5,932,177; 5,951,952; 6,267,929; 6,309,890 and 6,340,573; and Foreign Counterparts. Other Patents Pending.

BIOMÉRIEUX

The logo is a registered and protected trademark of bioMérieux sa or one of its subsidiaries / Le logo est une marque déposée et protégée par bioMérieux SA ou l'une de ses filiales / Das Logo ist ein eingetragenes und geschütztes Warenzeichen von bioMérieux sa bzw. deren Tochtergesellschaften / El logo es una marca registrada y protegida propiedad exclusiva de bioMérieux S.A. o de cada una sus filiales / Il logo é un marchio depositato e protetto di proprietà esclusiva di bioMérieux sa o di una delle sue filiali / O logótipo é uma marca comercial registada e legalmente protegida propriedade exclusiva da bioMérieux sa ou de uma das suas filiais / Dette logo er et registreret og beskyttet varemærke tilhørende bioMérieux sa eller et af dets datterselskaber / Logotypen är ett registrerat och skyddat varumärke som tillhör bioMérieux sa eller ett av dess dotterbolag / Το λογότυπο αποτελεί καταχωρημένο και προστατευόμενο εμπορικό σήμα της bioMérieux SA ή μίας εκ των θυγατρικών της / Logo jest zarejestrowanym i zastrzeżonym znakiem towarowym firmy bioMérieux sa lub jednego z jej przedstawicielstw

bioMérieux, Inc.
Box 15969, Durham NC 27704-0969 USA
tel. (1) 800 682 2666

EC|REP
bioMérieux® sa
au capital de 12 029 370 €
673 620 399 RCS LYON / 69280 Marcy-l'Etoile / France
tel. (33) 04 78 87 20 00 / fax (33) 04 78 87 20 90
http://www.biomerieux.com