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
T0 = cTn concentration on admission | T2h = cTn concentration 2 hours after admission

* Data from validation study in 592 patients | ** Difference in cTn concentration

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VIDAS® High sensitive Troponin I	
Reference number	415386
Tests / kit	60
Time to result	20 min
Sample type	Serum, Plasma (Li Hep)
Sample volume	200 µL
Units	ng/L
Measuring range	4.9 – 40 000 ng/L
Limit of quantification (LoQ)	2.9-4.9 ng/L
Limit of detection (LoD)	1.3-3.2 ng/L
Limit of blank (LoB)	0.0-1.9 ng/L
Precision (99 th percentile for the global healthy population)	19 ng/L ; CV = 7.0%
Calibrators & Controls frequency	28 days

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Did you know?

Worldwide **16 million** people die annually of **cardiovascular diseases**, specifically heart attacks or stroke ⁽¹⁾

Chest pain is among the **top 10** reasons to visit an emergency department (ED), accounting for about 5% of all ED visits ⁽²⁾


Expected prevalence of acute myocardial infarction (MI) in chest pain patients in the ED ⁽³⁾

5-10%
STEMI*

15-20%
NSTEMI**

Optimize the management of ACS*** with high sensitivity troponins ⁽⁴⁾

- New generation of **highly sensitive troponin (cTn) assays** allows **earlier detection of acute MI**, with shortening of time window for serial measurement to 3 hours.
- cTn has become a continuous variable with accurate measurement below the 99th percentile as well as small absolute changes within one or two hours. This has enabled development of **algorithms for reliable rule-out and rule-in of acute MI within 2 hours**.

Diagnosis

- Mandatory tool, in addition to clinical assessment and ECG, to make a definitive diagnosis of acute MI.
- Serial measurement to differentiate acute from chronic cardiac myocyte damage.

Risk stratification

Additional tool in risk stratification to guide therapeutic decision making.

Deliver high performance

Intended use:

- Aid in the **diagnosis of myocardial infarction**.
- Aid in the **risk stratification** of patients with symptoms suggestive of ACS with respect to relative risk of all-cause mortality and major adverse cardiac events (MACE) consisting of MI and revascularization, at 30 days.


The **99th** percentile has been defined at 19 ng/L

At the 99th percentile

Accelerate patient triage with high diagnostic accuracy

- Algorithm based on admission values and absolute changes within the first 2 hours.
- Delivers a diagnosis - safe rule-out and accurate rule-in - of acute MI (NSTEMI) in the vast majority (70%) of suspected chest pain patients.
- Shortens the time needed to triage chest pain patients in the ED (will obviate the need for prolonged monitoring with frequent serial blood sampling).

* ST - segment elevation myocardial infarction
** Non - ST - segment elevation myocardial infarction
*** Acute Coronary Syndrome

* Data from validation study in 592 patients
** Difference in cTn concentration

*** Negative Predictive Value
**** Predictive Positive Value