Intended use

For determination of resistance to activated Protein C, caused by the Factor V:Q506 (Factor V Leiden) mutation, in plasma from untreated individuals and from patients on oral anti-coagulant (OAT) or heparin therapy.

Summary and principle

The APC resistance phenotype^{1,2} is, in more than 90% of cases, due to a mutation in the Factor V gene, resulting in a replacement of Arg⁵⁰⁶ (R) with Gln (Q) in the Factor V protein.^{3,4,5} The selectivity for the Factor V:Q506 or other mutations in the Factor V gene rendering the protein resistant to inactivation by APC is increased by normalizing the concentrations of other plasma proteins involved in formation and regulation of thrombin. By performing the APTT-based APC resistance assay in the presence of an excess of Factor V Reagent Plasma, the sensitivity and specificity for the Factor V:Q506 mutation is significantly increased. Further, this modification allows for the analysis of plasma from patients who are on OAT.^{6,7,8,9,10,11} Sample plasma is prediluted with Factor V Reagent Plasma and incubated with the APTT reagent for a standard period of time. Coagulation is triggered by the addition of CaCl₂ in the absence and presence of APC and the time of clot formation is recorded.

Composition

The **Factor V Leiden (APC™ Resistance V)** kit consists of:

- R** **APTT reagent** (Cat. No. 0020008721): 2 x 4 mL vials of purified phospholipids with colloidal silica as contact activator.
- FV** **Factor V Reagent Plasma** (Cat. No. 0020008723): 2 x 4 mL vials of lyophilized human plasma with a low level of Factor V activity and filler (also available separately from Instrumentation Laboratory, Cat. No. 0020008800).
- A** **APC/Calcium chloride** (Cat. No. 0020008722): 2 x 2 mL vials of human activated Protein C co-lyophilized with CaCl₂.
- Ca** **Calcium chloride** (Cat. No. 0020008720): 2 x 2 mL vials of calcium chloride in Tris buffer containing bovine serum albumin.
- 1** **APC Control Plasma Level 1** (Cat. No. 0020008724): 2 x 1 mL vials of lyophilized human normal plasma.
- 2** **APC Control Plasma Level 2** (Cat. No. 0020008725): 2 x 1 mL vials of lyophilized human abnormal plasma.

***Note:** APTT reagent and APC/CaCl₂ are not interchangeable between lots.

PRECAUTIONS AND WARNINGS:

The material in this product was tested by FDA approved test methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV antibodies. Handle as if potentially infectious.¹² Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29). Wear suitable protective clothing (S 36).

This product is For *in vitro* Diagnostic Use.

Preparation

APTT reagent: Mix thoroughly before use.
Factor V Reagent Plasma: Dissolve the contents of each vial with 4 mL of NCCLS Type II water or equivalent.¹³ Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use. Do not shake.
APC/Calcium chloride: Dissolve the contents of each vial with 2 mL of NCCLS Type II water or equivalent.¹³ Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use. Do not shake.
Calcium chloride: Supplied ready for use.
APC Control Plasma Level 1 & 2: Dissolve the contents of vial with 1 mL of NCCLS Type II water or equivalent.¹³ Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the control at 15-25°C for 30 minutes and invert to mix before use. Do not shake.

Reagent storage and stability

Unopened reagents and controls are stable until the expiration date shown on the vial when stored at 2-8°C.
APTT reagent - Opened reagent is stable: 1 week at 15-25°C, 1 month at 2-8°C in the original vial, 8 hours at 15°C on the ACL[®] Systems or 3 days at 15°C on the ACL Futura[®]/ACL Advance Systems and ACL TOP[™]. No stirring is required. Do not freeze.
Factor V Reagent Plasma - Stability after reconstitution: 8 hours at 15-25°C, 24 hours at 2-8°C, 3 months at -20°C or below in the original vial or 3 days at 15°C on the ACL Futura/ACL Advance Systems and ACL TOP. No stirring is required. Frozen reagent should be thawed at 37°C and gently mixed before use. Do not refreeze.

APC/Calcium chloride - Stability after reconstitution: 8 hours at 15-25°C, 5 days at 2-8°C, 3 months at -20°C or below in the original vial or 3 days at 15°C on the ACL Futura/ACL Advance Systems and ACL TOP. No stirring is required. Do not refreeze.
Calcium chloride - Opened reagent is stable: 1 week at 15-25°C, 1 month at 2-8°C in the original vial or 3 days at 15°C on the ACL Futura/ACL Advance Systems and ACL TOP. No stirring is required.
APC Control Plasma Level 1 & 2 - Stability after reconstitution: 6 hours at 2-8°C and 15-25°C or 3 months at -20°C or below in the original vial. Frozen control should be thawed at 37°C and gently mixed before use. Do not refreeze.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vials.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to NCCLS Document H21-A4 for further instructions on specimen collection, handling and storage.¹⁴

Additional reagents and control plasma

The following may be purchased separately.

| | |
|-------------------------|----------------------------------|
| | Americas, Pacific Rim and Europe |
| Factor V Reagent Plasma | Cat. No. 0020008800 |
| Cleaning solution | Cat. No. 0009831700 |

Quality control

Normal and abnormal controls are recommended for a complete quality control program.¹⁵ APC Control Plasmas Level 1 (normal) and Level 2 (abnormal) are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.¹⁶

Results

APCR-V results may be reported in seconds and ratios. Refer to the instrument's Operator's Manual for additional information.

Limitations/interfering substances

No significant differences are obtained between fresh and frozen samples. The sensitivity and specificity for the Factor V:Q506 mutation on analysis of plasma from OAT patients is not affected by the INR value. No interference up to 1 U/mL (unfractionated and low molecular weight) heparin on the IL Coagulation Systems. APCR-V results on the ACL Futura/ACL Advance Systems are not affected by hemoglobin up to 400 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 1000 mg/dL. APCR-V results on the ACL TOP are not affected by hemoglobin up to 500 mg/dL, bilirubin up to 18.8 mg/dL and triglycerides up to 1791 mg/dL. Patients with known high inhibitor activity (e.g. phospholipid antibodies) may give an abnormal APTT and thus possibly misleading results. In such cases, increasing the dilution factor (e.g. 1+9 or 1+19) may correct the test result.

Cut-off value

APC Resistance due to FV:Q506 mutation is indicated when the APCR-V ratio is below or equal to the cut-off value. When using the Factor V Leiden (APC Resistance V) kit, a cut-off value must be established for each individual instrument. The cut-off value is that value which distinguishes the abnormal from normal samples. The abnormal samples have APCR-V ratios which are below or equal to the cut-off value. The APCR-V ratio for normals is higher than the cut-off value.

Procedure to establish the cut-off value:

- Using a plasma sample with a normal APCR-V genotype, perform five independent determinations of the APCR-V ratio. Perform testing in triplicate. Confirm that the inter and intra assay variation is below 7%.
- Determine the APCR-V ratio for at least 30 plasma samples from individuals in the age range of 20 - 65 years. Include APC Control Plasma Level 1 and Level 2 for assay validation.
- Verify that APCR-V ratios for the APC Control Plasmas are within the range on the insert sheet.
- Calculate the median APCR-V ratio of the 30 samples.


- Calculate the Factor V related APC Resistance cut-off value as follows: if the median is < 2.8: 0.8 x median APCR-V ratio (calculated in #4); if the median is > 2.8: 0.75 x median APCR-V ratio (calculated in #4).
- The APCR-V ratio for APC Control Plasma Level 1 should be higher than the established cut-off value. The APCR-V ratio for APC Control Plasma Level 2 should be lower than the calculated cut-off value. IL recommends to perform regular checks of the cut-off value by use of the APC Control Plasma Level 1 and Level 2.

Expected values

A normal range study was performed using Factor V Leiden (APC Resistance V) kit.

| System | N | Range (units) |
|------------------------|-----|-----------------------------|
| ACL Family | 30 | 2.18 - 3.38 (APCR-V ratios) |
| ACL Futura/ACL Advance | 30 | 2.42 - 3.18 (APCR-V ratios) |
| ACL TOP | 120 | 2.61 - 3.32 (APCR-V ratios) |

The following results illustrate the applicability of the APC Resistance V assay on the analysis of plasma from untreated individuals, patients on OAT with an INR of 1.3 - 6.0, and patients receiving unfractionated Heparin or low molecular weight heparin (shown here are Fragmin[®] and Fraxiparin[®]). By prediluting patient plasma with Factor V Reagent Plasma, a complete discrimination between normal, heterozygous V:Q506 and homozygous V:Q506 genotypes can be obtained by performing APC resistance analysis on plasma from both untreated patients and patients on oral anticoagulant or heparin therapy.

Sample APCR-V ratios are distributed according to the Factor V genotype. For a given genotype, similar APCR-V ratios are obtained regardless of patient group. Typical median values are 2.6 (normal), 1.7 (heterozygous FV:Q506) and 1.2 (homozygous FV:Q506). Factor V Leiden (APC Resistance V) kit provides > 99% sensitivity for FV:Q506 mutation as determined on the ACL system (n = 295). These results were obtained using a specific lot of reagent. Due to many variables which may affect clotting times, each laboratory should establish its own normal range.

Performance characteristics

Precision: Within run and total (run to run and day to day) precision were assessed over multiple runs using both normal and abnormal samples.

| ACL Family | Mean (APCR-V ratios) | CV % (Within run) | CV % (Total) |
|-------------------------------|-----------------------------|--------------------------|---------------------|
| APC Control 1 | 2.49 | 1.8 | 2.0 |
| APC Control 2 | 1.58 | 1.9 | 2.2 |
| ACL Futura/ACL Advance | Mean (APCR-V ratios) | CV % (Within run) | CV % (Total) |
| APC Control 1 | 2.34 | 2.1 | 2.3 |
| APC Control 2 | 1.54 | 1.2 | 1.1 |
| ACL TOP | Mean (APCR-V ratios) | CV % (Within run) | CV % (Total) |
| APC Control 1 | 2.94 | 3.9 | 4.6 |
| APC Control 2 | 1.85 | 1.4 | 2.0 |

| Correlation: Unit | system | slope | intercept | r | Reference method |
|-------------------|---------|-------|-----------|-------|-----------------------|
| APCR-V ratios | ACL TOP | 1.041 | -0.004 | 0.970 | APCR-V on ACL Advance |

These precision and correlation results were obtained using specific lots of reagents and controls.

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