Factor XIII Antigen - 0020201300

Intended use

Automated latex enhanced immunoassay for the quantitative determination of Factor XIII Antigen (FXIII Ag) in human citrated plasma on IL Coagulation Systems.

Summary and principle

Factor XIII, also called fibrin-stabilizing factor, is a zymogen composed of two potentially active A subunits and two B subunits (carrier protein) that circulates in plasma as a tetramer. In the terminal phase of the clotting cascade, by the action of thrombin and Ca²⁺, the A subunits dissociate from the B subunits assuming an enzymatically active configuration.¹⁻³ FXIII is essential for maintaining hemostasis due to its role in fibrin stabilization and in the protection of fibrin from proteolytic degradation by the fibrinolytic system.²

Patients with inherited FXIII deficiencies exhibit severe bleeding diathesis and, in most cases, require life-long supplementation therapy. In addition, FXIII is also involved in maintaining pregnancy and in wound healing.³

Acquired FXIII deficiency may occur in several diseases, including inflammatory bowel diseases and acute leukemia. Increased plasma FXIII activity has been reported in patients with obliterative atherosclerosis and diabetic angiopathy, and in chronic leukemia patients with increased megakaryocytic activity.

Recently, a polymorphism of the FXIII A subunit, has been associated with a protective effect against occlusive vascular diseases.²⁻³

The FXIII Ag Latex Reagent is a suspension of uniform size polystyrene latex particles coated with rabbit polyclonal antibodies, highly specific for the A-subunit of FXIII. When a plasma containing the active A-subunit of FXIII is mixed with the Latex Reagent and the Buffer included in the kit, the coated latex particles agglutinate. The degree of agglutination is directly proportional to the concentration of FXIII Ag in the sample and is determined by measuring the decrease of transmitted light caused by the aggregates.⁴

Composition

The Factor XIII Antigen kit consists of:

- R** **Latex Reagent** (Cat. No. 0020201310): 2 vials x 2.5 mL of a suspension of polystyrene latex particles coated with a rabbit polyclonal antibody directed against the A-subunit of FXIII containing buffer, stabilizer and preservative.
- B** **Buffer** (Cat. No. 0020201320): 2 vials x 5 mL of TRIS buffer containing bovine serum albumin, stabilizers and preservative.
- Di** **Diluent** (Cat. No. 0020201330): 2 vials x 6 mL of TRIS buffer containing bovine serum albumin, stabilizers and preservative.

PRECAUTIONS AND WARNINGS:

Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29). Wear suitable protective clothing (S 36).

All reagents contain less than 0.1% sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures.

This product is For in vitro Diagnostic Use.

Preparation

Invert to mix all reagents before use.

Note: Avoid foam formation when mixing reagents. Bubbles on top of the liquids may interfere with the instruments liquid sensors.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

Opened reagents are stable 1 month at 2-8°C in the original vial or 4 days at 15°C on the ACL Futura®/ACL Advance Systems and ACL TOP™. Do not freeze.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

NOTE: The Factor XIII Antigen assay should be run in the batch mode with a cleaning cycle performed before and after running the Antigen assay on the IL Coagulation Systems to eliminate any possible reagent carryover.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to NCCLS Document H21-A4 for further instructions on specimen collection, handling and storage.⁵ Thaw frozen specimens at 37°C for at least 15 minutes and centrifuge plasma before testing. After thawing the assay must be performed within 2 hours.

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim Cat. No.	Europe Cat. No.
Calibration plasma	0020003700	0020003700
Normal Control	0020003120	0020003110
Special Test Control Level 2	0020010200	0020010200

Quality control

Normal and abnormal controls are recommended for a complete quality control program.⁶ Normal Control and Special Test Control Level 2 are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard et al for identification and resolution of out-of-control situations.⁷

Results

FXIII Ag results are reported in % of normality. Refer to the instrument's Operator's Manual for additional information. The assay results should be used with other information, including the clinical context, in forming a diagnosis.

Limitations/interfering substances

FXIII Ag results on ACL Futura/ACL Advance and ACL TOP Systems are not affected by hemoglobin up to 500 mg/dL, bilirubin up to 18 mg/dL, triglycerides up to 1280 mg/dL and Rheumatoid Factor up to 500 IU/mL.

Expected values

A normal range study was performed using the FXIII Ag kit.

System	N	Range (units)
ACL Futura/ACL Advance	121	79.6 - 150.5 (% FXIII Ag)
ACL TOP	121	75.2 - 154.8 (% FXIII Ag)

Ranges were calculated as recommended by NCCLS Document C28-A2.⁸ These results were obtained using a specific lot of reagent. Due to many variables which may affect results, each laboratory should establish its own FXIII Ag normal range.

Performance characteristics

Precision: Within run and total (run to run and day to day) precision was assessed over multiple runs.

ACL Futura ACL Advance	Mean (% FXIII Ag)	CV% (Within run)	CV% (Total)
Normal Control	77.0	2.9	3.9
Special Test Control Level 2	36.1	4.9	7.5

ACL TOP	Mean (% FXIII Ag)	CV% (Within run)	CV% (Total)
Normal Control	75.2	2.9	3.5
Special Test Control Level 2	37.7	4.5	5.5

Correlation: System	slope	intercept	r	Comparative method
ACL Futura ACL Advance	0.948	5.987	0.976	HemosIL FXIII Ag on ACL TOP
ACL TOP	1.017	-7.788	0.950	FXIII photometric assay on photo optical instrument

The precision and correlation results were obtained using specific lots of reagents and controls.

Detection limit:

System	
ACL Futura/ACL Advance	1.8 (% FXIII Ag)
ACL TOP	2.5 (% FXIII Ag)

Linearity: System	Units	Regular Test	Rerun test (High range)	Rerun test (Low range)
ACL Futura ACL Advance	(% FXIII Ag)	15 - 150	45 - 450	2 - 50
ACL TOP	(% FXIII Ag)	15 - 150	30 - 300	3.8 - 45

Without prozone effect up to 870% on the ACL Futura/ACL Advance and ACL TOP.

If the linear range (Regular test) is exceeded at the lower end (i.e. FXIII Ag results below 15%) then samples could be reassayed triggering the Rerun test (Low range).

Instruments with Auto Rerun capability (i.e. ACL TOP) perform this new assay and result correction automatically thereby expanding the test range down to 3.8%.

If the linear range (Regular test) is exceeded at the higher end (i.e. FXIII Ag results above 150%) then samples will be reassayed triggering the Auto Rerun capability and the result will be corrected automatically thereby expanding the test range up to 300% (ACL TOP) or up to 450% (ACL Futura/ACL Advance).

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