**Fibrinogen-C - 0020301100**

**Intended use**

For the quantitative determination of fibrinogen, based on the Clauss method, in human citrated plasma on IL Coagulation Systems.

**Summary and principle**

Several congenital abnormalities of fibrinogen result in impaired conversion of fibrinogen to fibrin during blood coagulation.<sup>1</sup> Fibrinogen is also a useful marker in the evaluation of several disease states including Disseminated Intravascular Coagulation, liver disease, inflammatory diseases and malignancies.<sup>2</sup> High levels of fibrinogen are associated with an increased risk for cardiovascular disease.<sup>3,4,5</sup> Increased levels are also found during pregnancy and oral contraceptive use, while reduced levels are found during thrombolytic therapy.<sup>2</sup>

The Fibrinogen-C kit uses an excess of thrombin to convert fibrinogen to fibrin in diluted plasma. At high thrombin and low fibrinogen concentration, the rate of reaction is a function of fibrinogen concentration.<sup>9</sup>

**Composition**

The Fibrinogen-C kit consists of:

**T Bovine thrombin** (Cat. No. 0020301110): 10 x 2 mL vials of lyophilized bovine thrombin (35 UNIH/mL) with bovine albumin, calcium chloride, buffer and stabilizers.

**PRECAUTIONS AND WARNINGS:**

All animal products should be treated as potentially infectious. Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29). Wear suitable protective clothing (S 36). This product is For *in vitro* Diagnostic Use.

**Preparation**

**Bovine thrombin:** Dissolve the contents of each vial with 2 mL of CLSI (formerly NCCLS) Type CLRW water or equivalent.<sup>7</sup> Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use. Do not shake.

**Reagent storage and stability**

Unopened reagent is stable until the expiration date shown on the vial when stored at 2-8°C.

**Bovine Thrombin** - Stability after reconstitution: 3 days at 2-8°C, 1 month at -20°C in the original vial, 8 hours at 15°C on the ACL® Systems and 3 days at 15°C on the ACL Futura®/ACL Advance and ACL TOP® Family.

**Instrument/test procedures**

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

**Specimen collection and preparation**

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to CLSI Document H21-A5 for further instructions on specimen collection, handling and storage.<sup>8</sup>

**Additional reagents and control plasmas**

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim Cat. No.	Europe Cat. No.
Calibration plasma	0020003700	0020003700
Normal Control	0020003120	0020003110
Low Abnormal Control	0020003220	0020003210
Low Fibrinogen Control	0020004200	0020004200
Factor diluent	0009757600	0009757600
Cleaning solution	0009831700	0009831700
Cleaning agent	0009832700	0009832700

**Quality control**

Normal and abnormal controls are recommended for a complete quality control program.<sup>9,10</sup> Normal Control, Low Abnormal Control and Low Fibrinogen Control are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.<sup>11</sup>

**Results**

Patients results may be reported in the following units:  
- mg/dL or g/L  
- seconds  
Refer to the instrument's Operator's Manual for additional information.

**Limitations/interfering substances**

Fibrinogen-C results on the ACL and ACL Futura/ACL Advance Systems are not affected by heparin up to 1 U/mL, hemoglobin up to 150 mg/dL, triglycerides up to 1022 mg/dL or bilirubin up to 19.6 mg/dL. Fibrinogen-C results on ACL TOP Family are not affected by heparin up to 1 U/mL, hemoglobin up to 375 mg/dL, triglycerides up to 750 mg/dL, bilirubin up to 21 mg/dL. Fibrinogen-C assay results may be affected by degradation products (fibrin or fibrinogen) in the plasma assayed.<sup>12</sup>

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**ENGLISH - Insert revision 11/2009**

**Expected values**

A normal range study was performed using Fibrinogen-C reagent.

System	N	Range (unit)
ACL Familie	35	220 - 496 (mg/dL)
ACL Futura/ACL Advance	35	203 - 472 (mg/dL)
ACL TOP Familie	119	238 - 498 (mg/dL)

These results were obtained using a specific lot of reagent. Due to many variables which may affect clotting times, each laboratory should establish its own normal range.

**Performance characteristics**

**Precision:**

Within run and total (run to run and day to day) precision was assessed over multiple runs using both normal and abnormal samples.

ACL Family	Mean (Fibrinogen mg/dL)	CV % (Within run)	CV % (Total)
Normal Control	334	1,31	3,69
Low Fibrinogen Control	96	2,83	5,94

ACL Futura/ACL Advance	Mean (Fibrinogen mg/dL)	CV % (Within run)	CV % (Total)
Normal Control	329	2,75	4,93
Low Fibrinogen Control	101	3,68	4,94

ACL TOP Family	Mean (Fibrinogen mg/dL)	CV % (Within run)	CV % (Total)
Normal Control	303	4,5	5,5
Low Fibrinogen Control	107	5,1	6,8

**Correlation:**

System	slope	intercept	r	Reference method
ACL Futura/ACL Advance	0,94	-8,73	0,978	Fibrinogen-C on ACL
ACL Familie	0,88	22,90	0,983	Clauss Fibrinogen
ACL TOP Familie	0,98	-8,90	0,996	Clauss Fibrinogen on ACL ELITE/ELITE PRO

These precision and correlation results were obtained using specific lots of reagents and controls.

**Linearity:**

System	ACL TOP Familie	ACL Futura/ACL Advance	ACL Familie
	35 - 1000 mg/dL	80 - 700 mg/dL	70 - 700 mg/dL

Due to many variables which may affect results, each laboratory should establish its own Linearity Range.

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