

Heparin - 0020009400

Intended use

Automated chromogenic assay for the quantitative determination of unfractionated heparin (UFH) and low molecular weight heparin (LMWH) activity in human citrated plasma on the IL Coagulation Systems.

Summary and principle

Heparin is the most frequently used antithrombotic drug. The biological activity of this sulphated glycosaminoglycan resides in its ability to accelerate (up to 2000-fold) the inhibitory effect of antithrombin on coagulation proteases. In recent years, it has been shown that LMWH, besides being as useful therapeutically as UFH, also has a longer half-life.

The Heparin kit is an assay based on a synthetic chromogenic substrate and on Factor Xa inactivation. Heparin levels in patient plasma are measured automatically on IL Coagulation Systems in two stages.

- Heparin is analyzed as a complex with antithrombin present in the sample. The concentration of this complex is dependent on the availability of antithrombin. In order to obtain a more constant concentration of antithrombin, purified human antithrombin is added to the test plasma.¹ Factor Xa is added in excess and is neutralized by heparin-antithrombin complex.
- Residual Factor Xa is quantified with a synthetic chromogenic substrate. The paranitroaniline released is monitored kinetically at 405 nm and is inversely proportional to the heparin level in the sample.²

Since different kinds of UF- and LMW- Heparins have their own specific anti-Factor Xa activity, the same kind of heparin as is used in the patient sample should also be used for calibrating the standard curve.^{3,4}

Composition

The Heparin kit consists of:

- S** **Chromogenic substrate** (Cat. No. 00020009410): 1 x 4 mL vial of the lyophilized chromogenic substrate S-2765, N- α -Z-D-Arg-Gly-Arg-pNA-2HCl (3 mg/vial) and bulking agent.
- E** **Factor Xa reagent** (Cat. No. 00020009420): 1 x 5 mL vial of a lyophilized preparation containing purified bovine Factor Xa (68 nkat/vial), Tris-buffer, EDTA, sodium chloride and bovine serum albumin.
- A** **Antithrombin** (Cat. No. 00020009430): 1 x 3 mL vial of a lyophilized preparation containing human antithrombin (3 IU/vial), Tris-buffer, EDTA, sodium chloride and bovine serum albumin.
- B** **Buffer** (Cat. No. 00020009440): 1 x 8 mL vial of a concentrated solution containing Tris-buffer, pH 8.4, EDTA, sodium chloride and detergent.

PRECAUTIONS AND WARNINGS:

The material in this product was tested with FDA cleared methods and found nonreactive for Hepatitis B surface Antigen (HBsAg), Anti-HCV and HIV antibodies. Handle as if potentially infectious.⁵

Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29). Wear suitable protective clothing (S 36).

All animal products should be treated as potentially infectious.

This product is For *in vitro* Diagnostic Use.

Preparation

Chromogenic substrate: Dissolve the vial contents with 4 mL of NCCLS Type II water or equivalent.⁶ Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the substrate at 15-25°C for 30 minutes and invert to mix before use.

Factor Xa reagent: Dissolve the vial contents with 5 mL of NCCLS Type II water or equivalent.⁶ Replace the stopper and swirl gently. Make sure of the complete reconstitution of product. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use.

Antithrombin: Dissolve the vial contents with 3 mL of NCCLS Type II water or equivalent.⁶ Replace the stopper and swirl gently. Make sure of the complete reconstitution of product. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use.

Buffer: Dilute the necessary quantity of concentrated buffer 1:10 (1+9) with NCCLS Type II water or equivalent.⁶ Mix before use.

Working buffer: To 24 mL of diluted buffer add 1 mL of reconstituted antithrombin reagent.

Note: An instant opalescence will occur in the lyophilized reagents but it will fade away within 2 minutes.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

Chromogenic substrate - Stability after reconstitution: 7 days at 15°C, 3 months at 2-8°C in the original vial or 48 hours at 15°C on the ACL Futura®/ACL Advance Systems and ACL TOP™.

Factor Xa reagent - Stability after reconstitution: 7 days at 15°C, 3 months at 2-8°C in the original vial or 48 hours at 15°C on the ACL Futura®/ACL Advance Systems and ACL TOP.

Antithrombin - Stability after reconstitution: 3 months at 2-8°C in the original vial.

Buffer - Opened reagent is stable 3 months at 2-8°C.

Working buffer - Stability after preparation: 7 days at 15°C and 2-8°C in a closed container or 48 hours at 15°C on the ACL Futura®/ACL Advance Systems and ACL TOP.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

Note: A cleaning cycle is recommended after running the Heparin assay on all IL Coagulation Systems.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to NCCLS Document H21-A4 for further instructions on specimen collection, handling and storage.⁷

Standardization

For preparation of standards use the same heparin as used for patient therapy. Refer to the instrument's Operator's Manual for complete instructions on their preparation.

Additional reagent and control plasmas

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim Cat. No.	Europe Cat No.
Calibration plasma	0020000000	0008467300
Low Heparin control*	20004000	N/A
High Heparin control*	20004100	N/A
Cleaning solution	0009831700	0009831700
Cleaning agent	0009832700	0009832700

***Note:** Not for use with low molecular weight heparin (LMWH) assays.

Quality control

Two levels of controls are recommended for a complete quality control program.⁸ Low* and High* Heparin control are designed for this program.

Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution for out-of-control situations.⁹

Results

Heparin results are reported in U/mL. Refer to the instrument's Operator's Manual for additional information.

Limitations/interfering substances

Heparin results on the ACL and ACL Futura/ACL Advance Systems are not affected by hemoglobin up to 200 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 700 mg/dL.

Heparin results on the ACL TOP are not affected by hemoglobin up to 375 mg/dL, bilirubin up to 25 mg/dL and triglycerides up to 1630 mg/dL.

Expected values

To obtain an optimal effect with minimum risk of bleeding or thromboembolic complications the heparin activity should be in the range recommended by the heparin manufacturer.¹⁰

Performance characteristics

Precision:

Within run and total (run to run and day to day) precision was assessed over multiple runs using two levels of sample of both the 4th International Standard for heparin (UFH) and Fragmin® (LMWH).

ACL® Family	Mean (U/mL)	CV % (Within run)	CV % (Total)
UFH	0.77	1.84	2.18
UFH	0.23	7.76	8.23
LMWH	0.79	2.68	3.09
LMWH	0.23	7.99	9.69

ACL Futura/ACL Advance	Mean (U/mL)	CV % (Within run)	CV % (Total)
UFH	0.79	3.0	6.6
UFH	0.52	5.7	7.3
UFH	0.26	9.1	10.0
LMWH	0.76	4.1	4.5
LMWH	0.42	6.2	7.9
LMWH	0.22	10.7	11.9

ACL TOP	Mean (U/mL)	CV % (Within run)	CV % (Total)
UFH	0.82	1.8	6.4
UFH	0.53	3.4	4.3
UFH	0.25	4.6	8.5
LMWH	0.86	2.4	4.4
LMWH	0.49	6.1	8.2
LMWH	0.25	5.7	11.3

Correlation:

System	slope	intercept	r	Reference method
ACL Family	0.968	0.014	0.988	IL Heparin (Xa)
ACL Futura/ACL Advance	0.944	0.035	0.989	IL Heparin (Xa)
ACL TOP	1.012	-0.005	0.992	HemosIL Heparin (Xa) on ACL Advance

The precision and correlation results were obtained using specific lots of reagent and controls.

Linearity:

System






ACL Family and ACL Futura/ACL Advance	0 - 1.0 U/mL
ACL TOP	0 - 1.1 U/mL

Bibliography / Literatur / Bibliografia / Bibliographie / Bibliografía / Bibliografia / Litteratur / Litteraturförteckning / Βιβλιογραφία

- Teien AN, Lie M. Evaluation of an Amidolytic Heparin Assay Method: Increased Sensitivity by Adding Purified Antithrombin III, *Thromb. Res.* 1977; 10: 399-410.
- Teien AN, *et al.* Assay of Heparin in Plasma using a Chromogenic Substrate, *Thromb. Res.* 1976; 8: 413-416.
- Fareed J, Hoppenstadt D, Jeske W, Clarizio R, Waleng JM, Bick RL. The Available Low Molecular Weight Heparin Preparations Are Not the Same, *Clin. Appl. Thrombosis/Hemostasis.* 1997; 3 (suppl. 1): s38-s52.
- Holmer E, Kurachi K, Söderström. The Molecular-Weight Dependence of the Rate-Enhancing Effect of Heparin on the Inhibition of Thrombin, Factor Xa, Factor IXa and Kallikrein by Antithrombin, *Biochem J.* 1981: 193: 395-400.
- Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.
- National Committee for Clinical Laboratory Standards. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, NCCLS Document C3-A3; Vol. 17 No. 18.
- National Committee for Clinical Laboratory Standards. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, NCCLS Document H21-A4; Vol. 23 No. 35.
- Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, *Am J. Clin. Pathol.* 1970; 53: 924-927.
- Westgard JO, and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
- Holm HA, Abildgaard U, Kalvenes S. Heparin Assays and Bleeding Complications in Deep Venous Thrombosis with Particular Reference to Retroperitoneal Bleeding, *Thromb. Haemost.* 1985; 53: 278-281.

ACL, ACL Futura and ACL TOP are trademarks of Instrumentation Laboratory
 ©1998 Instrumentation Laboratory
 Issued March 2005

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<p>IVD</p> <p><i>In vitro</i> diagnostic medical device <i>In-vitro</i> Diagnostikum De uso diagnóstico <i>in vitro</i> Dispositif médical de diagnostic <i>in vitro</i> Per uso diagnostico <i>in vitro</i> Dispositivo médico para utilização em diagnóstico <i>in vitro</i> "in vitro" diagnostisk udstyr <i>In vitro</i> diagnostisk medicinsk produkt</p> <p>Προϊόν για διαγνωστική χρήση <i>In vitro</i></p>	<p>LOT</p> <p>Batch code Chargen-Bezeichnung Identificación número de lote Désignation du lot Numero del lotto Número de lote Batch nr. Tillverkningskod Αρ. Παρτίδας</p>	<p></p> <p>Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzare prima del Data limite de utilização Anvendelse Användning Χρήση έως</p>	<p></p> <p>Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Limite de temperatura Temperatur begrænsninger Temperatur gräns Περιορισμοί θερμοκρασίας</p>	<p></p> <p>Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização Se vejledning for anvendelse Ta del av instruktionen före användning Συμβουλευτήτε τις οδηγίες χρήσης</p>	<p>CONTROL</p> <p>Control Kontrollen Control Contrôle Controllo Controllo Kontrol Kontroll Υλικό ποιστικού ελέγχου</p>	<p></p> <p>Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico Miljø oplysninger Biologiska risker Βιολογικοί κίνδυνοι</p>	<p></p> <p>Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por Producent Tillverkare Κατασκευαστής</p>	<p>EC REP</p> <p>Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado Leverandør Auktoriserad representant Εξουσιοδοτημένος αντιπρόσωπος</p>
--	--	--	---	---	---	---	--	---