# Liquid Heparin - 0020300100

### Intended use

Automated chromogenic assay for the quantitative determination of unfractionated heparin (UFH) and low molecular weight heparin (LMWH) activity in human citrated plasma on IL Coagulation Systems, (ACL TOP™ Family, ACL ELITE /ELITE PRO/8/9/10000 and ACL Futura™/ACL Advance™).

### Summary and principle

Heparin is the most frequently used antithrombotic drug. The biological activity of this sulphated glycosaminoglycan resides in its ability to accelerate (up to 2000-fold) the inhibitory effect of antithrombin on coagulation proteases. In recent years, it has been shown that LMWH, besides being as useful therapeutically as UFH, also has a longer half-life.

The Liquid Heparin kit is a one stage chromogenic assay based on a synthetic chromogenic substrate and on Factor Xa inactivation. Heparin levels in patient plasma are measured automatically on IL Coagulation Systems.

Heparin is analyzed as a complex with antithrombin present in the sample. The concentration of this complex is dependent on the availability of the patient's endogenous antithrombin. When the Heparin – antithrombin complex is formed, two competing reactions take place.

- Factor Xa is neutralized by heparin-antithrombin complex.
- Residual Factor Xa is quantified with a synthetic chromogenic substrate. The paranitroaniline released is monitored kinetically at 405 nm and is inversely proportional to the heparin level in the sample.<sup>1</sup>

In order to reduce the influence from heparin antagonists, such as platelet factor 4 (PF4), dextran sulfate is included in the reaction mixture.

### Composition

The **Heparin** kit consists of:

**[S] Chromogenic substrate** (Cat. No. 0020300120): 5 x 3 mL vial of liquid chromogenic substrate S-2732 (approximately 1.2 mg/mL) and bulking agent.

**[E] Factor Xa reagent** (Cat. No. 0020300110): 5 x 2.5 mL vial of a liquid preparation containing purified bovine Factor Xa (approximately 5.5 nkat/mL), Tris-Buffer, EDTA, dextran sulfate, sodium chloride and bovine serum albumin.

#### PRECAUTIONS AND WARNINGS:

Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29). Wear suitable protective clothing (S 36).

All animal products should be treated as potentially infectious.

This product is For *in vitro* Diagnostic Use.

#### Preparation

**Chromogenic substrate:** Invert to mix before use.

**Factor Xa reagent:** Invert to mix before use.

#### Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

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**Chromogenic substrate** – Opened reagent is stable: 1 month at 2-8°C or 7 days at 15° -25°C on-board the ACL TOP Family, 3 days on-board the ACL ELITE/ELITE PRO/ 8/9/10000 and 4 days on-board the ACL Futura/ACL Advance Systems in the original vial.

**Factor Xa reagent** - Opened reagent is stable: 1 month at 2-8°C or 7 days at 15° -25°C on-board the ACL TOP Family, 3 days on-board the ACL ELITE/ELITE PRO/8/9/10000 and 4 days on-board the ACL Futura/ACL Advance Systems in the original vial.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

### Instrument/test procedures

Refer to the appropriate IL instrument’s Operator’s Manual for the complete assay procedure instructions.

### Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to CLSI (formerly NCCLS) Document H21-A5 for further instructions on specimen collection, handling and storage.<sup>2</sup>

### Additional reagent and control plasmas

The following are not supplied with the kit and must be purchased separately.

Heparin Calibrators	Product	Cat No.
LMW Heparin Controls	0020300600	
UF Heparin Controls	0020300200	
Cleaning solution	0009831700	
Cleaning agent	0009832700	

### Quality control

Two levels of controls are recommended for a complete quality control program.<sup>3</sup> LMW and UF Heparin Controls are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument’s Operator’s Manual for additional information. Refer to Westgard *et al* for identification and resolution for out-of-control situations.<sup>4</sup>

### Results

Heparin results are reported in IU/mL. Refer to the instrument’s Operator’s Manual for additional information.

### Limitations/interfering substances

Heparin results on the ACL TOP Family, ACL ELITE/ELITE PRO/8/9/10000 and ACL Futura/ACL Advance Systems are not affected by hemoglobin up to 300 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 800 mg/dL.

### Expected values

To obtain an optimal effect with minimum risk of bleeding or thromboembolic complications the heparin activity should be in the range recommended by the heparin manufacturer.<sup>5</sup>

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### Performance characteristics

#### Precision

Within run and total (run - to - run and day - to - day) precision was assessed over multiple runs using the two levels of both UFH and LMWH Controls.

ACL TOP Family	Mean (IU/mL)	CV% (Within run)	CV% (Total)
UFH Low	0.41	2.4	3.3
UFH High	0.68	1.0	1.6
LMWH Low	0.55	3.5	4.5
LMWH High	1.35	1.9	2.5
ACLELITE/ELITE PRO/8/9/10000	Mean (IU/mL)	CV% (Within run)	CV% (Total)
UFH Low	0.39	4.0	6.7
UFH High	0.69	1.0	3.7
LMWH Low	0.49	5.5	5.7
LMWH High	1.31	3.3	4.0
ACL Futura/ACL Advance	Mean (IU/mL)	CV% (Within run)	CV% (Total)
UFH Low	0.41	4.4	4.4
UFH High	0.69	1.3	1.7
LMWH Low	0.56	2.5	3.5
LMWH High	1.37	1.1	1.6

<b>Correlation:</b>				
System	Slope	Intercept	r	Reference method
ACL TOP Family	0.946	0.055	0.958	HemosIL Heparin
ACL ELITE/ELITE PRO/8/9/10000	0.894	0.057	0.907	HemosIL Heparin
ACL Futura/ACL Advance	1.067	0.006	0.946	HemosIL Heparin

The precision and correlation results were obtained using specific lots of reagent and controls.

In field site studies using samples from patients undergoing heparin therapy, the following data were obtained using a specific lot of Liquid Heparin reagents:

System	Slope	Intercept	r	Reference method
ACL TOP Family	0.952	0.095	0.978	HemosIL Heparin
ACL ELITE/ELITE PRO/8/9/10000	1.032	0.039	0.949	HemosIL Heparin
ACL Futura/ACL Advance	1.007	-0.002	0.957	HemosIL Heparin

#### Detection limit:

System	
ACL TOP Family	0.04 IU/mL
ACL ELITE/ELITE PRO/8/9/10000	0.04 IU/mL
ACL Futura/ACL Advance	0.02 IU/mL

#### Linearity:

System	
ACL TOP Family	Up to - 2.0 IU/mL
ACL ELITE/ELITE PRO/8/9/10000	Up to - 2.0 IU/mL
ACL Futura/ACL Advance	Up to - 2.0 IU/mL

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