Intended use

Diluted Russell's Viper Venom Test (DRVVT) reagents for the detection of lupus anticoagulants (a type of phospholipid interfering antibody) in human citrated plasma on the IL Coagulation Systems.
LAC Screen: Simplified DRVVT reagent to screen for the presence of Lupus Anticoagulants.
LAC Confirm: Phospholipid-rich DRVVT reagent to confirm the presence of Lupus Anticoagulants.
 These tests are indicated for use with patients who have a prolonged APTT test of undetermined origin.

Summary and principle

Lupus anticoagulants (LA) belong to the group of antiphospholipid antibodies which are directed against negatively charged phospholipids or complexes between phospholipids and proteins (either beta-2-glycoprotein 1 or clotting factors such as prothrombin). When determined by their ability to prolong phospholipid-dependent tests (APTT, KCT, DRVVT) these antibodies are referred to as LA. Patients with LA are at increased risk of clinical complications such as thrombosis and recurrent miscarriages.^{1,2,3,4}
 LAC Screen and LAC Confirm are improved DRVVT reagents and intended to simplify and standardize the detection of LA in clinical evaluations.¹ LAC Screen is poor in phospholipid making it sensitive to LA. The additional amount of phospholipid in LAC Confirm neutralizes LA to give shorter clotting times. Russell's viper venom, in the presence of calcium, directly activates factor X (in the test sample). LAC Screen and LAC Confirm are therefore unaffected by contact factor abnormalities, factor VII, VIII and IX deficiencies or inhibitors. Heparin interference up to 1 U/mL is neutralized by polybrene. As a result, LAC Screen and LAC Confirm are more specific tests for the evaluation of LA than APTT.⁵

Composition

Each LAC kit consists of:

- LS** **LAC Screen** (Cat. No. 0020008000): 10 x 2 mL vials of a lyophilized preparation containing Russell's viper venom, phospholipids, calcium, polybrene, buffers, stabilizers, dyes and preservative.
- LC** **LAC Confirm** (Cat. No. 0020008200): 10 x 2 mL vials of a lyophilized preparation containing Russell's viper venom, phospholipids, calcium, polybrene, buffers, stabilizers, dyes and preservative.

PRECAUTIONS AND WARNINGS:

All reagents contain less than 0.1% sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures.
 Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29).
 Wear suitable protective clothing (S 36).

This product is For in vitro Diagnostic Use.

Preparation

Dissolve the contents of each vial with 2 mL of NCCLS Type II water or equivalent.⁶ Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use. Do not shake.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.
 Stability after reconstitution: 24 hours at 15-25°C, 48 hours at 2-8°C, 1 month at -20°C in the original vial, 2 days at 15°C on the ACL® Futura/ACL Advance Systems or 3 days at 15°C on the ACL TOP™. No stirring is required. Do not refreeze.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Specimen collection, preparation and handling

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to NCCLS Document H21-A4 for further instructions on specimen collection, handling and storage.⁷
 Freezing of the plasma releases residual platelet phospholipids which can shorten LAC Screen clotting times. Double centrifuge or filter plasma through 0.2 micron filter to remove platelets before freezing.

Instrument/test procedures

It is recommended that the LAC Screen and LAC Confirm be used together.

ACL® (100 - 7000) Method

1. From the PROG screen select Acquisition Time. Select Thrombin Time and enter into extended time.
2. Load the sample tray with sample plasmas.
3. Add LAC Screen or LAC Confirm reagent to reagent position 1.
4. From The Main Test Menu screen, select Thrombin Time and start analysis.
5. Sample results will be displayed and printed out (in seconds).
6. Sample results must be used in a calculation for the appropriate ratio result to be reported.

ACL Futura/ACL Advance, ACL 8000/9000/10000 and ACL TOP Method

Refer to the appropriate IL instrument's Operator's Manual for the complete assay procedure instructions.

Quality control

Normal and abnormal controls are recommended for a complete quality control program.⁸ Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard et al for identification and resolution of out-of-control situations.⁹

Results

The following procedure should be used to calculate the LAC Screen and LAC Confirm ratios:

1. For each new lot of LAC Screen and LAC Confirm kit a new Normal Range should be determined according to NCCLS Document H21-A4, Vol. 23, No. 35.
2. Determine the Mean of each Normal Range in seconds.
3. The mean of each normal range will be used as a constant denominator in the calculations of ratios.

LAC Screen

1. The patient sample result in seconds is divided by the Mean of the LAC Screen normal range.

$$\text{Screen Ratio} = \frac{\text{Patient Screen results (in seconds)}}{\text{Mean of Screen Normal Range (in seconds)}}$$

2. If the LAC Screen clotting time is within the normal range no further testing for LA may be necessary.
3. If test plasma LAC Screen clotting time is 20% longer than the Mean of the Screen Normal Range (ie ratio >1.2), the presence of LA should be confirmed with LAC Confirm.

LAC Confirm

1. The patient sample result in seconds is divided by the Mean of the LAC Confirm normal range.

$$\text{Confirm Ratio} = \frac{\text{Patient Confirm results (in seconds)}}{\text{Mean of Confirm Normal Range (in seconds)}}$$

2. The ratio results from the LAC Screen is divided by the ratio result from LAC Confirm.

$$\text{Normalized LAC Ratio} = \frac{\text{Screen Ratio}}{\text{Confirm Ratio}}$$

Interpretation

1. The final result should be expressed as the Normalized LAC Ratio:
 Ratio greater than 2.0 LA is strongly present.
 Ratio between 1.5 - 2.0 LA is moderately present.
 Ratio between 1.2 - 1.5 LA is weakly present.
2. If ratio < 1.2 and LAC Screen and LAC Confirm clotting times are prolonged, then mixing studies should be performed. Mixing normal plasma with test plasma replenishes any factors deficient in the test plasma.

Limitations/interfering substances

Icteric, lipemic and hemolyzed specimens may give false results.
 LAC Screen/LAC Confirm results on the ACL Futura/ACL Advance Systems are not affected by hemoglobin up to 400 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 295 mg/dL.
 LAC Screen/LAC Confirm results on the ACL Systems are not affected by hemoglobin up to 100 mg/dL, bilirubin up to 21.7 mg/dL and triglycerides up to 329 mg/dL.
 LAC Screen/LAC Confirm results on the ACL TOP are not affected by hemoglobin up to 125 mg/dL and bilirubin up to 13 mg/dL. Lipemic or turbid samples should not be assayed.
 LA assays based on different properties appear to be more or less sensitive to certain subgroups of LAs. Therefore at least two screening assays, based on different properties, should be performed before the possibility of LA is excluded.¹⁰

Expected values

A normal range study was performed using LAC Screen/LAC Confirm on the IL Coagulation Systems.

System	Normalized LAC Ratio
ACL Family and ACL Futura/ACL Advance	0.8 - 1.2
ACL TOP	0.9 - 1.1

These results were obtained using a specific lot of reagent. Due to many variables which may affect clotting times, each laboratory should establish its own normal range.

Performance characteristics

Within run and total (run to run and day to day) precision was assessed over multiple runs using both normal and abnormal samples.

ACL Family	Mean (Normalized LAC Ratio)	CV % (Within run)	CV % (Total)
Normal sample	0.93	1.16	1.70
Weakly LA sample	1.33	0.84	3.02
Strongly LA sample	2.03	3.63	4.55

ACL Futura/ACL Advance	Mean (Normalized LAC Ratio)	CV % (Within run)	CV % (Total)
Normal sample	1.00	1.1	2.2
Weakly LA sample	1.31	1.2	1.9
Strongly LA sample	1.81	1.7	2.2

ACL TOP	Mean (Normalized LAC Ratio)	CV % (Within run)	CV % (Total)
Normal sample	0.97	1.6	2.5
Weakly LA sample	1.30	1.9	2.4
Strongly LA sample	1.82	2.7	2.9

Normalized LAC ratio > 1.2 was found in the following plasma samples using LAC Screen/LAC Confirm:

known LA	90% (26/29)
heparin	12% (1/8)
warfarin	0% (0/7)
normal	2% (1/60)
factor deficiencies	0% (0/5)

Correlation:

Unit	system	slope	intercept	r	Reference method
Normalized LAC ratio	ACL Futura/ACL Advance	0.922	0.048	0.988	LAC Screen/Confirm on ACL Classic
Normalized LAC ratio	ACL 8000/9000/10000	0.926	0.075	0.994	LAC Screen/Confirm on ACL Classic
Normalized LAC ratio	ACL TOP	1.009	-0.022	0.996	LAC Screen/Confirm on ACL Advance

The precision and correlation results were obtained using specific lots of reagents.

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1. Love PE, Santoro SA. Antiphospholipid Antibodies: Anti Cardiolipin and the Lupus Anticoagulant in Systemic Lupus Erythematosus (SLE) and in Non-SLE Disorders, Ann. Int. Med. 1990; 112: 682-698.
2. Exner T, Triplett DA, Taberner D, Machin SJ. Guidelines for Testing and Revised Criteria for Lupus Anticoagulants, Thromb. Haemostas. 1991; 65: 320-322.
3. Rauch J, Tannebaum M and Janoff AS. Distinguishing Plasma Lupus Antifactor Antibodies using Hexagonal Phase II Phospholipids, Thromb. Haemostas. 1989; 62: 892-896.
4. Feinstein DI. Lupus Anticoagulant, Thrombosis and Fetal Loss, N. Eng. J. Med. 1985; 313: 1248-1350.
5. Thiagarajan P, Pengo V, Shapiro SS. The use of the Dilute Russell's Viper Venom Time for the Diagnosis of Lupus Anticoagulants, Blood Coag. Fibrinol. 1990; 1: 259-266.
6. National Committee for Clinical Laboratory Standards. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, NCCLS Document C3-A3; Vol. 17 No. 18.
7. National Committee for Clinical Laboratory Standards. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays, Fourth Edition, NCCLS Document H21-A4; Vol. 23 No. 35.
8. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimen for Coagulation, Am. J. Pathol. 1970; 53: 924-927.
9. Westgard JO and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
10. Brandt JT, Triplett DA, Alving B, Scharer I. Criteria for the Diagnosis of Lupus Anticoagulants: An Update. Thromb.Haemostas. 1995, 74 (4), 1185-90

ACL, ACL Futura and ACL TOP are trademarks of Instrumentation Laboratory Issued December 2008

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<p>IVD</p> <p>In vitro diagnostic medical device In-vitro Diagnostikum De uso diagnóstico in vitro Dispositif médical de diagnostic in vitro Per uso diagnostico in vitro Dispositivo médico para utilização em diagnóstico in vitro "in vitro" diagnostisk udstyr In vitro diagnostisk medicinsk produkt</p> <p>Προϊόν για διαγνωστική χρήση In vitro</p>	<p>LOT</p> <p>Batch code Chargen-Bezeichnung Identificación número de lote Désignation du lot Numero del lotto Número de lote Batch nr. Tillverkningskod Αρ. Παρτίδας</p>	<p></p> <p>Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzare prima del Data limite de utilização Anvendelse Användning Χρήση έως</p>	<p></p> <p>Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Limite de temperatura Temperatur begrænsninger Temperatur gräns Περιορισμοί θερμοκρασίας</p>	<p></p> <p>Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização Se vejledning for anvendelse Ta del av instruktionen före användning Συμβουλευτήτε τις οδηγίες χρήσης</p>	<p>CONTROL</p> <p>Control Kontrollen Control Contrôle Controllo Controllo Kontrol Kontroll Υλικό ποιοτικού ελέγχου</p>	<p></p> <p>Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico Miljø oplysninger Biologiska risker Βιολογικοί κίνδυνοι</p>	<p></p> <p>Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por Producent Tillverkare Κατασκευαστής</p>	<p>EC REP</p> <p>Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado Leverandör Auktoriserad representant Εξουσιοδοτημένος αντιπρόσωπος</p>
---	--	---	--	--	---	--	---	---

LAC Screen/Confirm 20008000/20008200

Printed Insert Sheet: 301233
Revision: R5
Issued: December 2008
C.O.: 408136

LANGUAGES

ENGLISH
GERMAN
SPANISH
FRENCH
ITALIAN
PORTUGUESE
DANSK
SVENSK
GREEK

TECHNICAL SPEC'S

PAPER: White paper, 50-60 g/m² weight.
SIZE: N° 4 Sheets 11 x 8.5" (280 x 216 mm.).
PRINT: Front/Back.
PRINT COLOR: Front - Top band Green Pantone 382, all remaining type in black.
Back - All type in black.