

**Intended use**

Automated chromogenic assay for the quantitative determination of Plasminogen in human citrated plasma on IL Coagulation Systems.

**Summary and principle**

Plasminogen (Plg) in a plasma sample is activated through reaction with an excess of Streptokinase (Sk) in the presence of fibrinogen. The Plasminogen-Streptokinase complex is determined by the rate of hydrolysis of a chromogenic substrate. Dysplasminogenemia is associated with venous thrombosis. Acquired deficiencies of Plasminogen are associated with thrombolytic therapy, sepsis and Disseminated Intravascular Coagulation.

The Plasminogen kit is based on a synthetic chromogenic substrate. Plasminogen levels in patient plasma are measured automatically on IL Coagulation Systems in two stages:

1. Incubation of the plasma with Streptokinase reagent in the presence of fibrinogen.
2. Quantification of the Plasminogen-Streptokinase complexes with a synthetic chromogenic substrate. The paranitroaniline released is monitored kinetically at 405 nm and is directly proportional to the Plasminogen level in the test sample.

**Composition**

The Plasminogen kit consists of:

- S Chromogenic substrate** (Cat. No. 0020009010): 2 x 2 mL vials of the lyophilized chromogenic substrate S-2403, pyroGlu-Phe-Lys-pNAHCl (4.2 mg/vial) and bulking agent.
- E Streptokinase reagent** (Cat. No. 0020009020): 2 x 2.5 mL vials of a lyophilized preparation containing Streptokinase (20,000 U/vial), fibrinogen, buffer and human serum albumin.

**PRECAUTIONS AND WARNINGS:**

The material in this product was tested with FDA cleared methods and found nonreactive for Hepatitis B surface Antigen (HBsAg), Anti-HCV and HIV antibodies. Handle as if potentially infectious.

Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29). Wear suitable protective clothing (S 36).

This product is For *in vitro* Diagnostic Use.

**Preparation**

**Chromogenic substrate:** Dissolve the contents of each vial with 2 mL of NCCLS Type II water or equivalent. Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the substrate at 15-25°C for 30 minutes and mix gently before use.

**Streptokinase reagent:** Dissolve the contents of each vial with 2.5 mL of NCCLS Type II water or equivalent. Replace the stopper and swirl gently. Make sure of the complete reconstitution of the product. Keep the reagent at 15-25°C for 30 minutes and mix gently before use.

**Reagent storage and stability**

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

**Chromogenic Substrate** - Stability after reconstitution: 7 days at 15°C, 3 months at 2-8°C, 6 months at -20°C in the original vial or 5 days at 15°C on the ACL Futura /ACL Advance Systems.

**Streptokinase reagent** - Stability after reconstitution: 7 days at 15°C, 3 months at 2-8°C, 6 months at -20°C in the original vial or 5 days at 15°C on the ACL Futura/ACL Advance Systems.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

**Instrument/test procedures**

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

**Specimen collection and preparation**

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to NCCLS Document H21-A3 for further instructions on specimen collection, handling and storage.

**Additional reagents and control plasmas**

The following are not supplied with the kit and must be purchased separately.

	Americas and Pacific Rim Cat. No.	Europe Cat No.
Calibration plasma	0020000000	0008467300
Normal Control	0020003120	0020003110
Abnormal chromogenic control plasma Level 1/2	0008467600	0008467600
Factor Diluent	0009757600	0009757600

**Quality control**

Normal and abnormal controls are recommended for a complete quality control program. Normal Control and Abnormal chromogenic control plasma Level 1/2 are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.

**Results**

Plasminogen results are reported in % activity. Refer to the instrument's Operator's Manual for additional information.

**Limitations/interfering substances**

Plasminogen results on IL Coagulation Systems are not affected by elevated levels of fibrinogen, heparin (UF or LMW heparin) up to 2 U/mL, hemoglobin up to 200 mg/dL, bilirubin up to 20 mg/dL, FDP's up to 30 mg/dL and triglycerides up to 1000 mg/dL.

**Expected values**

A normal range study was performed using Plasminogen kit.

System	N	Range (units)
ACL Family	48	72.9 - 126.9 (% activity)
ACL Futura/ACL Advance	59	81.4 - 140.9 (% activity)

Ranges were calculated as recommended by the International Federation of Clinical Chemistry (IFCC). These results were obtained using a specific lot of reagent. Due to many variables which may affect results, each laboratory should establish its own normal range.

**Performance characteristics**

**Precision:**

Within run and total (run to run and day to day) precision was assessed over multiple runs using both normal and abnormal samples.

ACL Family	Mean (% activity)	CV % (Within run)	CV % (Total)
Normal Control	86.6	2.0	3.8
Abnormal Control	50.1	4.4	4.4

ACL Futura/ ACL Advance	Mean (% activity)	CV % (Within run)	CV % (Total)
Normal Control	82.7	3.0	3.0
Abnormal Control	52.5	1.8	1.9

Correlation: System	slope	intercept	r	Reference method
ACL Family	1.218	-6.064	0.983	Chromogenic Plasminogen
ACL Futura/ ACL Advance	1.119	-2.818	0.989	Chromogenic Plasminogen

The precision and correlation results were obtained using specific lots of reagents and controls.

**Linearity:  
System**

ACL Family and ACL Futura/ACL Advance 10 - 150 (% activity)

**Tests on the IL Coagulation Systems**

Plasminogen60 tests (approx.)






## Bibliography / Literatur / Bibliografía / Bibliographie / Bibliografia / Bibliografia / Litteratur / Litteraturförteckning / Βιβλιογραφία

1. Aoki N, Moroi M, Yoshida N, Matsuda M. Abnormal Plasminogen. A Hereditary Molecular Abnormality Found in a Patient with Recurrent Thrombosis, J. Clin. Invest. 1978; 61: 1186-1195.
2. Robbins KC. Dysplasminogenemias, Prog. Cardiovasc. Disease. 1992; 4: 295-308.
3. Latallo ZS, Teisseyre E, Lopaciuk S. Evaluation of a Fibrinolytic Profile of Plasma using Chromogenic Substrates. In: Chromogenic Peptide Substrate. Chemistry and Clinical usage. Scully MK, Kakkar VV (eds), Churchill Livingstone, London. 1979; 262-268.
4. Gallimore MJ, Aasen AO, Erichsen NS, Larsbraaten M, Lyngaas K, Amundsen E. Plasminogen Concentrations and Functional Activities and Concentration of Plasmin Inhibitors in Plasma Samples from Normal Subjects and Patients with Septic Shock, Thrombosis Research. 1980; 18: 601-608.
5. Soria J, Soria C, Samama M. A Plasminogen Assay using a Chromogenic Synthetic Substrate. Results from Clinical Work and from Studies of Thrombolysis. In: Progress of Chemical Fibrinolysis and Thrombolysis. Davidson JF, Rowan RM, Samama MM, Desnoyers PC (eds). Raven Press, New York. 1978; 337-346.
6. Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4<sup>th</sup> Edition, 1999.
7. National Committee for Clinical Laboratory Standards. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, NCCLS Document C3-A3; Vol. 17 No. 18.
8. National Committee for Clinical Laboratory Standards. Collection, Transport and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays, Third Edition, NCCLS Document H21-A3; Vol. 18 No. 20.
9. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am J. Clin. Pathol. 1970; 53: 924-927.
10. Westgard JO, and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
11. IFCC Committee Report, J. Clin. Chem. Clin Biochem. 1983; 21: 749-760.

ACL and ACL FUTURA are registered trademarks of Instrumentation Laboratory

©1998 Instrumentation Laboratory  
Issued July 2003

## Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<p><b>IVD</b></p> <p><i>In vitro</i> diagnostic medical device <i>In-vitro</i> Diagnostikum De uso diagnóstico <i>in vitro</i> Dispositif médical de diagnostic <i>in vitro</i> Per uso diagnostico <i>in vitro</i> Dispositivo médico para utilização em diagnóstico <i>in vitro</i> "in vitro" diagnostisk udstyr <i>In vitro</i> diagnostisk medicinsk produkt Προϊόν για διαγνωστική χρήση <i>In vitro</i></p>	<p><b>LOT</b></p> <p>Batch code Chargen-Bezeichnung Identificación número de lote Désignation du lot Numero del lotto Número de lote Batch nr. Tillverkningskod Αρ. Παρτίδας</p>	<p></p> <p>Use by Verwendbar bis Caducidad Utilisable jusqu'à Da utilizzzare prima del Data límite de utilização Anvendelse Användning Χρήση έως</p>	<p></p> <p>Temperature limitation Festgelegte Temperatur Temperatura de Almacenamiento Températures limites de conservation Limiti di temperatura Límite de temperatura Temperatur begrænsninger Temperatur gräns Περιορισμοί θερμοκρασίας</p>	<p></p> <p>Consult instructions for use Beilage beachten Consultar la metódica Lire le mode d'emploi Vedere istruzioni per l'uso Consultar as instruções de utilização Se vejledning for anvendelse Ta del av instruktionen före användning Συμβουλευτήτε τις οδηγίες χρήσης</p>	<p><b>CONTROL</b></p> <p>Control Kontrollen Control Contrôle Controllo Controlo Kontrol Kontroll Υλικό ποιοτικού ελέγχου</p>	<p></p> <p>Biological risks Biologisches Risiko Riesgo biológico Risque biologique Rischio biologico Risco biológico Miljø oplysninger Biologiska risker Βιολογικοί κίνδυνοι</p>	<p></p> <p>Manufacturer Hergestellt von Fabricado por Fabricant Prodotto da Fabricado por Producent Tillverkare Κατασκευαστής</p>	<p><b>EC REP</b></p> <p>Authorised representative Bevollmächtigter Representante autorizado Mandataire Rappresentanza autorizzata Representante autorizado Leverandør Auktoriserad representant Εξουσιοδοτημένος αντιπρόσωπος</p>
--	--	---	---	---	--	---	--	---