Silica Clotting Time - 0020004800

Intended Use

For the detection of Lupus Anticoagulants in human citrated plasma on the IL Coagulation System by the use of screening (SCT Screen) and confirmatory (SCT Confirm) reagents sensitized to phospholipid dependent antibodies.

Summary and Principle

Lupus anticoagulants (LA) belong to the group of antiphospholipid antibodies which are directed against negatively charged phospholipids or complexes between phospholipids and proteins (for example, beta-2-glycoprotein I or clotting factors such as prothrombin). They determine by their ability to prolong phospholipid-dependent tests (APTT, KCT, DRVVT, SCT) these antibodies are referred to as LA. Patients with LA are at increased risk of clinical complications such as thrombosis and recurrent miscarriages.¹⁻⁴ SCT Screen and SCT Confirm are reagents intended to simplify and standardize the detection of LA in clinical evaluations.⁵⁻⁷ SCT Screen is poor in phospholipid making it sensitive to LA. The addition of an amount of phospholipid in SCT Confirm neutralizes LA to give shorter clotting times.

Silica Clotting Time in the presence of calcium, directly activates the intrinsic pathway of coagulation. SCT Screen and SCT Confirm are therefore unaffected by factor VII deficiencies or inhibitors. Heparin interference up to 0.4 U/mL is neutralized by polybrene. Using a ratio of screen and confirm allows the SCT to be insensitive to warfarin treated samples.⁸ As a result, SCT Screen and SCT Confirm are more specific tests for the evaluation of LA than APTT or dilute PT.⁹

Composition

The Silica Clotting Time kit consists of:

- Sc** **SCT Screen** (Cat. No. 0020004821): 3 x 5 mL vials of a liquid preparation containing colloidal silica, buffer and preservative.
- Cf** **SCT Confirm** (Cat. No. 0020004822): 3 x 5 mL vials of a liquid preparation containing colloidal silica, synthetic phospholipids, buffer and preservative.
- Ca** **SCT CaCl₂** (Cat. No. 0020004823): 3 x 10 mL vials of Calcium Chloride (0.025 Mol/L) with polybrene and preservative.

PRECAUTIONS AND WARNINGS: Avoid contact with skin and eyes (S 24/25). Do not empty into drains (S 29). Wear suitable protective clothing (S 36). Calcium Chloride contains sodium azide that may form explosive azides in metal plumbing. Use proper disposal procedures. This product is For *in vitro* Diagnostic Use.

Preparation

SCT Confirm: Shake Silica dispersion vigorously for approximately 15 seconds or vortex for 5 seconds before use. Subsequent handling after initial mixing requires only gentle inversion 3-4 times.

SCT Screen: Add 50 µL of mixed SCT Confirm reagent to one vial of SCT Screen reagent. Shake Silica dispersion vigorously for approximately 15 seconds or vortex for 5 seconds before initial mixing. Subsequent handling after initial mixing requires only gentle inversion 3-4 times.

SCT CaCl₂: invert to mix before use.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

SCT Screen and Confirm: Opened reagent is stable 20 days at 2-8°C in the original vial or 5 days at 15°C on IL Coagulation System.

No stirring is required. Do not freeze.

SCT CaCl₂: Opened reagent is stable 20 days at 2-8°C in the original vial.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual and/or Application Manual.

Specimen collection, preparation and handling

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to CLSI (formerly NCCLS) Document H21-A4 for further instructions on specimen collection, handling and storage.¹⁰

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- Love PE, Santoro SA. Antiphospholipid Antibodies: Anti Cardiolipin and the Lupus Anticoagulant in Systemic Lupus Erythematosus (SLE) and in Non-SLE Disorders. *Int. Med.* 1992; 112: 682-698.
- Exner T, Triplett DA, Taberner D, Machin SJ. Guidelines for Testing and Revised Criteria for Lupus Anticoagulants. *Thromb. Haemostas.* 1991; 65: 320-322.
- Rauch J, Tannebaum M and Janoff AS. Distinguishing Plasma Lupus Antifactor Antibodies Using Hexagonal Phase II Phospholipids. *Thromb. Haemostas.* 1989; 62: 692-696.
- Feinstein DL. Lupus Anticoagulant, Thrombosis and Fetal Loss. *N. Eng. J. Med.* 1985; 313: 1248-1350.
- Chantarangkul V, Tripodi A, Arbini A, Mannucci PM. Silica clotting time (SCT) as a screening and confirmatory test for detection of the lupus anticoagulant. *Thromb. Res.* 1992; 67: 351-354.
- Chantarangkul V, Tripodi A, Clerici M, Bressi C, Mannucci PM. Laboratory Diagnosis of Lupus Anticoagulants: Effect of Residual Platelets in Plasma. Assessed by Staclot LA and Silica Clotting Time. *Thromb. Haemost.* 2002; 87: 854-8.

Freezing of the plasma releases residual platelet phospholipids which can shorten SCT Screen clotting times. Double centrifuge or filter plasma through 0.2 micron filter to remove platelets before freezing.

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately.

Americas and Pacific Rim	Europe
Cat. No. 0020003120	Cat. No. 0020003110

Quality control

Normal and Lupus controls are recommended for a complete quality control program.¹¹ Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hours shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.¹²

Results

The following procedure should be used to calculate the SCT Screen and SCT Confirm ratios:

- For each individual instrument and for each specific lot of SCT Screen and SCT Confirm kit, determine a new Normal Range according to CLSI Document C28-A2, Vol. 20, No. 13.¹³
- Determine the Mean of each Normal Range in seconds.
- The mean of each normal range will be used as a constant denominator in the calculations of ratios.

SCT Screen

- The patient sample result in seconds is divided by the Mean of the SCT Screen normal range.

$$\text{Screen Ratio} = \frac{\text{Patient Screen results (in seconds)}}{\text{Mean of Screen Normal Range (in seconds)}}$$

SCT Confirm

- The patient sample result in seconds is divided by the Mean of the SCT Confirm normal range.

$$\text{Confirm Ratio} = \frac{\text{Patient Confirm results (in seconds)}}{\text{Mean of Confirm Normal Range (in seconds)}}$$

- The ratio result from the SCT Screen is divided by the ratio result from SCT Confirm.

$$\text{Normalized SCT Ratio} = \frac{\text{Screen Ratio}}{\text{Confirm Ratio}}$$

Interpretation

- The final result should be expressed as the Normalized SCT Ratio:

System

ACL TOP	Ratio > 1.16	LA is present
ACL Futura/ACL Advance	Ratio > 1.24	LA is present
ACL ELITE/ELITE PRO 8/9/10000	Ratio > 1.20	LA is present

- If ratio ≤ 1.16 (ACL TOP), ≤ 1.24 (ACL Futura/ACL Advance) or ≤ 1.20 (ACL ELITE/ ELITE PRO/8/9/10000) and SCT Screen and SCT Confirm clotting times are prolonged, then mixing studies should be performed to investigate factor deficiencies or inhibitors. If the mixing test is still prolonged, it indicates that some anticoagulant other than LA may be present in the test plasma.

Limitations/interfering substances

SCT Screen/SCT Confirm results on the ACL Futura/ACL Advance and ACL ELITE/ELITE PRO 8/9/10000 are not affected by bilirubin up to 30 mg/dL, triglycerides up to 500 mg/dL and Heparin up to 0.4 U/mL. Do not use hemolyzed samples.

SCT Screen/SCT Confirm results on the ACL TOP are not affected by bilirubin up to 30 mg/dL, triglycerides up to 850 mg/dL, UF Heparin up to 0.5 U/mL and LMW Heparin up to 1.0 U/mL. Do not use hemolyzed samples.

LA assays based on different properties appear to be more or less sensitive to certain subgroups of LAs. Therefore at least two screening assays, based on different properties, should be performed before the possibility of LA is excluded.¹⁴

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