



ARCHITECT

SYSTEM

en

Free T₃

REF 7K63

48-9164/R3

B7K630

Read Highlighted Changes
Revised December, 2009


Free T₃



Customer Service: Contact your local representative or find country specific contact information on www.abbottdiagnostics.com

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Key to symbols used

REF	List Number	CONTROL NO.	Control Number
IVD	<i>In Vitro</i> Diagnostic Medical Device	REAGENT LOT	Reagent Lot
LOT	Lot Number	SAMPLE CUPS	Sample Cups
	Expiration Date	SEPTUM	Septum
	Store at 2-8°C	REACTION VESSELS	Reaction Vessels
	Consult instructions for use	REPLACEMENT CAPS	Replacement Caps
	Manufacturer	SN	Serial Number
		WARNING: SENSITIZER	Warning: May cause an allergic reaction

See **REAGENTS** section for a full explanation of symbols used in reagent component naming.

NAME

ARCHITECT Free T₃

INTENDED USE

The ARCHITECT Free T₃ (FT₃) assay is a Chemiluminescent Microparticle Immunoassay (CMIA) for the quantitative determination of free triiodothyronine (Free T₃) in human serum and plasma.

SUMMARY AND EXPLANATION OF TEST

3,5,3' Triiodothyronine (T₃) is a thyroid hormone with a molecular weight of 651 daltons¹ and a half-life in serum of 1.5 days.² T₃ circulates in the blood as an equilibrium mixture of free and protein bound hormone.³ T₃ is bound to thyroxine binding globulin (TBG), prealbumin, and albumin. The actual distribution of T₃ among these binding proteins is controversial as estimates range from 38-80% for TBG, 9-27% for prealbumin, and 11-35% for albumin.⁴ The binding of these proteins is such that only 0.2-0.4% of the total T₃ is present in solution as unbound or free T₃.⁵ This free fraction represents the physiologically active thyroid hormone.³

Free T₃ is typically elevated to a greater degree than free thyroxine (T₄) in Graves' disease.^{6,7} Occasionally, free T₃ alone is elevated (T₃ thyrotoxicosis) in about 5% of the hyperthyroid population.⁸ In contrast, levels of free T₄ are elevated to a greater degree than free T₃ in toxic multinodular goiter and excessive T₄ therapy.⁹ Serum free T₃ is useful in distinguishing these forms of hyperthyroidism. Free T₃ may also be important in monitoring patients on anti-thyroid therapy where treatment is focused on reducing the T₃ production and the T₄ conversion to T₃. Serum free T₃ may also be useful in assessing the severity of the thyrotoxic state.

The ARCHITECT Free T₃ assay is to be used as an aid in the assessment of thyroid status.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT Free T₃ assay is a two-step immunoassay to determine the presence of free (unbound) T₃ in human serum and plasma using Chemiluminescent Microparticle Immunoassay (CMIA) technology with flexible assay protocols, referred to as Chemiflex.

In the first step, sample and anti-T₃ coated paramagnetic microparticles are combined. Free T₃ (unbound) present in the sample binds to the anti-T₃ coated microparticles. After washing, T₃ acridinium labeled conjugate is added in the second step. Pre-Trigger and Trigger Solutions are then added to the reaction mixture; the resulting chemiluminescent reaction is measured as relative light units (RLUs). An inverse relationship exists between the amount of Free T₃ in the sample and the RLUs detected by the ARCHITECT *i* optical system.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Reagent Kit, 100 Tests/500 Tests

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT *i* Systems. Please contact your local distributor.

ARCHITECT Free T₃ Reagent Kit (7K63)

- **MICROPARTICLES** 1 or 4 Bottle(s) (6.6 mL/27.0 mL) anti-T₃ (sheep) coated Microparticles in MES buffer with sheep IgG stabilizers. Preservative: antimicrobial agent.
- **CONJUGATE** 1 or 4 Bottle(s) (5.9 mL/26.3 mL) T₃ acridinium-labeled Conjugate in citrate buffer with NaCl and Triton X-100 stabilizers. Minimum concentration: 0.33 ng/mL. Preservative: antimicrobial agent.

Other Reagents

ARCHITECT *i* Pre-Trigger Solution

- **PRE-TRIGGER SOLUTION** Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

ARCHITECT *i* Trigger Solution

- **TRIGGER SOLUTION** Trigger Solution containing 0.35N sodium hydroxide.

ARCHITECT *i* Wash Buffer

NOTE: Bottle and volume vary based on order.

- **WASH BUFFER** Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

WARNINGS AND PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use.
- Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Safety Precautions

- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹⁰ Biosafety Level 2¹¹ or other appropriate biosafety practices^{12,13} should be used for materials that contain or are suspected of containing infectious agents.

The following warnings and precautions apply to these components:

- **Microparticles**
- **Conjugate**



WARNING: H317 Contains methylisothiazolones. May cause an allergic skin reaction.

Prevention

P261 Avoid breathing mist / vapours / spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves / protective clothing / eye protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P333+P313 If skin irritation or rash occurs: Get medical advice / attention.
P363 Wash contaminated clothing before use.

This material and its container must be disposed of in a safe way.

- For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Handling Precautions

- Do not use reagent kits beyond the expiration date.
- **Do not mix reagents from different reagent kits.**
- Prior to loading the ARCHITECT Free T₃ Reagent Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**
- To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
- Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
- Over time, residual liquids may dry on the septum surface. These are typically dried salts which have no effect on assay efficacy.
- For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Storage Instructions

- 
- The ARCHITECT Free T₃ Reagent Kit must be stored at 2-8°C and may be used immediately after removal from 2-8°C storage.
 - When stored and handled as directed, reagents are stable until the expiration date.
 - The ARCHITECT Free T₃ Reagent Kit may be stored on-board the ARCHITECT *i* System for a maximum of 30 days. After 30 days, the reagent kit must be discarded. For information on tracking on-board time, refer to the ARCHITECT System Operations Manual, Section 5.
 - Reagents may be stored on or off the ARCHITECT *i* System. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** After reagents are removed from the system, you must initiate a scan to update the on-board stability timer.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results may be invalid and may require retesting. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

- The ARCHITECT Free T₃ assay file must be installed on the ARCHITECT *i* System from the ARCHITECT *i* Assay CD-ROM prior to performing the assay. For detailed instructions on assay file installation and on viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.
- For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.
- For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.
- The default result unit for the ARCHITECT Free T₃ assay is pg/mL. An alternate result unit, pmol/L, may be selected for reporting results by editing assay parameter "Result concentration units" to pmol/L. The conversion factor used by the system is 1.536.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

- Human serum (including serum collected in serum separator tubes) or plasma collected in sodium heparin, lithium heparin, or potassium EDTA anticoagulant tubes may be used in the ARCHITECT Free T₃ assay. Other anticoagulants have not been validated for use with the ARCHITECT Free T₃ assay. Follow the manufacturer's processing instructions for serum or plasma collection tubes.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.
- The ARCHITECT *i* System does not provide the capability to verify specimen type. It is the responsibility of the operator to verify the correct specimen types are used in the ARCHITECT Free T₃ assay.
- Use caution when handling patient specimens to prevent cross contamination. Use of disposable pipettes or pipette tips is recommended.
- Do not use heat-inactivated specimens.
- For optimal results, inspect all samples for bubbles. Remove bubbles with an applicator stick prior to analysis. Use a new applicator stick for each sample to prevent cross contamination.
- For optimal results, serum and plasma specimens should be free of fibrin, red blood cells or other particulate matter.
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.

- If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator or red blood cells. Specimens may be stored for up to 6 days at 2-8°C prior to being tested. If testing will be delayed more than 6 days, specimens should be frozen at -10°C or colder. Specimens stored frozen at -10°C or colder for 6 days showed no performance difference.
- Multiple freeze-thaw cycles of specimens should be avoided. Specimens must be mixed THOROUGHLY after thawing, by LOW speed vortexing or by gently inverting, and centrifuged prior to use to remove red blood cells or particulate matter to ensure consistency in the results.
- When shipped, specimens must be packaged and labeled in compliance with applicable state, federal and international regulations covering the transport of clinical specimens and infectious substances. Prior to shipment, it is recommended that specimens be removed from the clot, serum separator or red blood cells.

PROCEDURE

Materials Provided

- 7K63 ARCHITECT Free T₃ Reagent Kit

Materials Required but not Provided

- ARCHITECT *i* System
- ARCHITECT *i* Assay CD-ROM
- 7K63-01 ARCHITECT Free T₃ Calibrators
- 7K63-10 ARCHITECT Free T₃ Controls
- ARCHITECT *i* **PRE-TRIGGER SOLUTION**
- ARCHITECT *i* **TRIGGER SOLUTION**
- ARCHITECT *i* **WASH BUFFER**
- ARCHITECT *i* **REACTION VESSELS**
- ARCHITECT *i* **SAMPLE CUPS**
- ARCHITECT *i* **SEPTUM**
- ARCHITECT *i* **REPLACEMENT CAPS**
- For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

Assay Procedure

- Before loading the ARCHITECT Free T₃ Reagent Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that have settled during shipment:
 - Invert the microparticle bottle 30 times.
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
 - Once the microparticles have been resuspended, remove and discard the cap. Wearing clean gloves, remove a septum from the bag. Carefully snap the septum onto the top of the bottle.
 - **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
- Order tests.
- Load the ARCHITECT Free T₃ Reagent Kit on the ARCHITECT *i* System. Verify that all necessary assay reagents are present. Ensure that septums are present on all reagent bottles.
- The minimum sample cup volume is calculated by the system and is printed on the Orderlist report. No more than 10 replicates may be sampled from the same sample cup. To minimize the effects of evaporation verify adequate sample cup volume is present prior to running the test.
 - Priority: 75 µL for the first Free T₃ test plus 25 µL for each additional Free T₃ test from the same sample cup
 - ≤ 3 hours on-board: 150 µL for the first Free T₃ test plus 25 µL for each additional Free T₃ test from the same sample cup
 - > 3 hours on-board: additional sample volume is required. Refer to the ARCHITECT System Operations Manual, Section 5 for information on sample evaporation and volumes.

- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- ARCHITECT Free T₃ Calibrators and Controls should be mixed by gentle inversion prior to use.
- To obtain the recommended 150 µL volume requirements for the ARCHITECT Free T₃ Calibrators and Controls, hold the bottles **vertically** and dispense 4 drops of each calibrator or 4 drops of each control into each respective sample cup.
- Load samples.
- For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN. The ARCHITECT *i* System performs the following functions:
 - Moves the sample to the aspiration point
 - Loads a reaction vessel (RV) into the process path
 - Aspirates and transfers sample into the RV
 - Advances the RV one position and transfers microparticles into the RV
 - Mixes, incubates and washes the reaction mixture
 - Adds conjugate to the RV
 - Mixes, incubates and washes the reaction mixture
 - Adds Pre-Trigger and Trigger Solutions
 - Measures chemiluminescent emission to determine the quantity of Free T₃ in the sample
 - Aspirates contents of RV to liquid waste and unloads RV to solid waste
 - Calculates the result
- For information on ordering patient specimens, calibrators and controls, and general operating procedures refer to the ARCHITECT System Operations Manual, Section 5.
- For optimal performance, it is important to follow the routine maintenance procedures defined in the ARCHITECT System Operations Manual, Section 9. If your laboratory requires more frequent maintenance, follow those procedures.

Specimen Dilution Procedures

Specimens cannot be diluted for Free T₃ determinations. Specimens which read > 30.00 pg/mL should be reported as such.

Calibration

- To perform an ARCHITECT Free T₃ calibration, test Calibrators 1 and 2 in duplicate. A single sample of all levels of Free T₃ controls must be tested to evaluate the assay calibration. Ensure that assay control values are within the concentration ranges specified in the package insert. Calibrators should be priority loaded.
- Calibrator Range: 0.0 - 30.0 pg/mL.
- Once an ARCHITECT Free T₃ calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used
 - Controls are out of range
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

QUALITY CONTROL PROCEDURES

The recommended control requirement for the ARCHITECT Free T₃ assay is a single sample of all control levels tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures. Ensure that assay control values are within the concentration ranges specified in the package insert.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B. The ARCHITECT Free T₃ assay belongs to method group 2.

RESULTS

The ARCHITECT Free T₃ assay utilizes a 4 Parameter Logistic Curve Fit data reduction method (4PLC, Y weighted) to generate a calibration curve.

Alternate Result Units

- The default result unit for the ARCHITECT Free T₃ assay is pg/mL. When the alternate result unit, pmol/L, is selected, the conversion factor used by the system is 1.536.
- Conversion Formula:
(Concentration in pg/mL) x (1.536) = Concentration in pmol/L

Flags

- Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other thyroid tests, clinical impressions, etc.
- If the Free T₃ results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- Performance of this test has not been established with neonatal specimens.

EXPECTED VALUES

A normal range of 1.71 pg/mL to 3.71 pg/mL (central 95% interval) was obtained by testing serum specimens from 436 individuals determined as normal by AxSYM Ultrasensitive hTSH II and AxSYM Free T₄ assays. It is recommended that each laboratory establish its own normal range, which may be unique to the population it serves depending upon geographical, patient, dietary, or environmental factors.

Free T₃ is a secondary indicator of thyroid status. Although the majority of patients with hyperthyroidism will have free T₃ values greater than the upper limit of the euthyroid range, some may have free T₃ values which fall within the normal range.^{14,15} Specimens from patients described as “sick euthyroids” generally yield values in the low to normal range.^{16,17}

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The ARCHITECT Free T₃ assay is designed to have a precision of ≤ 10% (total CV). A study based on guidance from Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) document EP5-A¹⁸ was performed for the ARCHITECT Free T₃ assay. A three member processed human serum based panel was assayed, using two lots of reagents, in replicates of two at two separate times per day for 20 testing days. Data from this study are summarized in the following table.*

Panel Member	Reagent Lot	Instrument	n	Mean Conc. Value (pg/mL)	SD	Within Run %CV	Total SD	Total %CV
1	1	1	80	3.22	0.096	3.0	0.115	3.6
1	1	2	80	3.16	0.133	4.2	0.143	4.5
1	2	1	80	3.60	0.108	3.0	0.141	3.9
1	2	2	80	3.35	0.113	3.4	0.131	3.9
2	1	1	80	6.00	0.099	1.7	0.168	2.8
2	1	2	80	5.88	0.166	2.8	0.184	3.1
2	2	1	80	6.28	0.154	2.5	0.176	2.8
2	2	2	80	6.06	0.194	3.2	0.225	3.7
3	1	1	80	10.50	0.252	2.4	0.481	4.6
3	1	2	80	10.01	0.289	2.9	0.496	5.0
3	2	1	80	10.50	0.145	1.4	0.237	2.3
3	2	2	80	10.12	0.217	2.1	0.265	2.6

* Representative data; results in individual laboratories may vary from these data.

Analytical Sensitivity

The ARCHITECT Free T₃ assay is designed to have an analytical sensitivity of ≤ 1.0 pg/mL. Analytical sensitivity is defined as the concentration calculated as the mean plus two standard deviations of replicates of the ARCHITECT Free T₃ MasterCheck Level 0 (0.0 pg/mL). The analytical sensitivity (low-linearity) is defined in the ARCHITECT Free T₃ assay parameters as 1.0 pg/mL.

Analytical Specificity

The ARCHITECT Free T₃ assay is designed to have a mean analytical specificity of $\leq 0.001\%$ cross reactivity with thyroxine (T₄) at a concentration of 1,000,000 pg/mL.

Interference

The ARCHITECT Free T₃ assay is designed to have a mean potential interference from hemoglobin, bilirubin, triglycerides, and protein of $< 10\%$ at the levels indicated below.

- Hemoglobin - ≤ 500 mg/dL
- Bilirubin - ≤ 20 mg/dL
- Triglycerides - ≤ 2000 mg/dL
- Protein - ≤ 12 g/dL

Accuracy by Correlation

The ARCHITECT Free T₃ assay is designed to have a slope of 1.00 +/- 0.20 and a correlation coefficient (r) of ≥ 0.90 when compared to the AxSYM Free T₃ assay. A study was performed where specimens were tested using the ARCHITECT Free T₃ assay and AxSYM Free T₃ assay. Data from this study were analyzed using least squares and Passing Bablok¹⁹ regression methods and are summarized in the following table.*

Abbott ARCHITECT Free T₃ vs. Abbott AxSYM Free T₃

Method	Number of Specimens	Intercept	Slope	Correlation Coefficient
Least Squares Linear Regression	1387	0.53	0.96	0.957
Passing-Bablok Linear Regression**	1387	0.24	1.07	0.957

In this evaluation, serum specimens tested ranged from 1.10 pg/mL to 29.91 pg/mL with the ARCHITECT Free T₃ assay and from 1.10 pg/mL to 27.57 pg/mL with the AxSYM Free T₃ assay.

* Representative data; variables such as differences in sampling size and sample population may impact the correlation of the assay; therefore, results in individual laboratories may vary from these data.

** A linear regression method with no special assumptions regarding the distribution of the samples and the measurement errors.¹⁹

BIBLIOGRAPHY

1. Budavari S, editor. *Merck Index* (11th Ed.). Rahway, NJ: Merck and Co., Inc., 1989:868.
2. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972; 21:1073-1092.
3. Ekins RP, editor. *Methods for the Measurement of Free Thyroid Hormones*. Amsterdam: Excerpta Medica Foundation. 1979;72-92.
4. Robbins J, Rall JE. The Iodine-Containing Hormones. In: *Hormones in Blood* (3rd Ed.). London: Academic Press, 1979:1:632-667.
5. DeGroot LJ, Larsen PR, Refetoff S, Stanbury JB. Transport of Thyroid Hormone and Cell Uptake. In: *The Thyroid and Its Diseases*. New York: Wiley and Sons, 1984:62-66.
6. Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously Functioning Thyroid Nodules. *J Clin Endocrinol Metab* 1980;50: 1089-1093.
7. Ladenson PW. Diagnosis of Thyrotoxicosis. In: Braverman LE, Utiger RD, editors. *The Thyroid* (6th Ed.). Philadelphia: JB Lippincott Co., 1991:880-886.
8. Wahner HW. T₃ Hyperthyroidism. *Mayo Clin Proc* 1972;47:938-943.
9. Lum SM, Nicoloff JT, Spencer CA, Kaptein EM. Peripheral Tissue Mechanism for Maintenance of Serum Triiodothyronine Values in a Thyroxine-Deficient State in Man. *J Clin Invest* 1984;73:570-575.
10. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.

11. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; January 2007.
12. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
13. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
14. Lee LA, Mooney RA, Woolf PD. Clinical Utility of Measuring Free Thyroxin and Free Triiodothyronine in Serum of Critically Ill Patients by Ultrafiltration. *Clin Chem* 1986; 32:797-800.
15. Wilke TJ, Eastment HT. Discriminative Ability of Tests for Free and Total Thyroid Hormones in Diagnosing Thyroid Disease. *Clin Chem* 1986; 32:1746:50.
16. Kaptein EM, Grieb DA, Spencer CA, et al. Thyroxine Metabolism in the Low Thyroxine State of Critical Nonthyroidal Illnesses. *J Clin Endocrinol Metab* 1981; 53:764-71.
17. Merimee TJ, Fineberg ES. Starvation-Induced Alterations of Circulating Thyroid Hormone Concentrations in Man. *Metabolism* 1976; 25:79-83.
18. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices - Approved Guideline. NCCLS Document EP5-A. Wayne, PA: NCCLS, 1999.
19. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. *J Clin Chem. Clin Biochem.* 1983;21:709-20.

ARCHITECT, AxSYM, Chemiflex and MasterCheck are trademarks of Abbott Laboratories in various jurisdictions.

All trademarks are property of their respective owners.



Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000



Distributed by Abbott Laboratories
Abbott Park, IL 60064 USA
and
ABBOTT 65205 Wiesbaden, Germany

December 2009
© 2005, 2009 Abbott Laboratories