



# ARCHITECT

## SYSTEM

# en

Free T<sub>4</sub>

**REF** 7K65

48-9249/R3

**B7K650**

Read Highlighted Changes  
Revised January, 2010

# Free T<sub>4</sub>



Customer Service: Contact your local representative or find country specific contact information on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

### Key to symbols used

|   |   |                            |   |
|---|---|----------------------------|---|
| <b>REF</b>  | List Number                               | <b>CAL 1</b>               | Calibrator (1, 2)   |
| <b>IVD</b>  | <i>In Vitro</i> Diagnostic Medical Device | <b>SEPTUM</b>              | Septum  |
| <b>LOT</b>  | Lot Number                                | <b>SN</b>                  | Serial Number   |
|  | Expiration Date                           | <b>CONTROL L</b>           | Control Low, Medium, High (L, M, H)                                 |
|  | Store at 2-8°C                            | <b>SAMPLE CUPS</b>         | Sample Cups   |
|  | Caution                                   | <b>REAGENT LOT</b>         | Reagent Lot   |
|  | Consult instructions for use              | <b>ASSAY CD-ROM</b>        | Assay CD-ROM  |
|  | Manufacturer                              | <b>REACTION VESSELS</b>    | Reaction Vessels  |
|   |   | <b>REPLACEMENT CAPS</b>    | Replacement Caps  |
|   |   | <b>WARNING: SENSITIZER</b> | Warning: May cause an allergic reaction                             |
|   |   | <b>CONTAINS: AZIDE</b>     | Contains sodium azide. Contact with acids liberates very toxic gas. |

See **REAGENTS** section for a full explanation of symbols used in reagent component naming.

## NAME

ARCHITECT Free T<sub>4</sub>

## INTENDED USE

The ARCHITECT Free T<sub>4</sub> (FT4) assay is a Chemiluminescent Microparticle Immunoassay (CMIA) for the quantitative determination of free thyroxine (Free T<sub>4</sub>) in human serum and plasma.

## SUMMARY AND EXPLANATION OF TEST

Thyroxine (T<sub>4</sub>) circulates in the blood as an equilibrium mixture of free and serum protein bound hormone. Thyroxine binding globulin (TBG), albumin and pre-albumin bind approximately 75%, 10% and 15% of the total circulating T<sub>4</sub> respectively.<sup>1-3</sup> The binding of T<sub>4</sub> by these proteins is such that less than 0.03% is present in the circulation as unbound, free T<sub>4</sub>.<sup>4</sup> This small percentage of the total T<sub>4</sub> represents the physiologically available hormone which is biologically active. Once the free T<sub>4</sub> is absorbed by the target cells, the equilibrium reestablishes circulating free T<sub>4</sub> levels. The equilibrium results in the maintenance of a constant level of free T<sub>4</sub> when alterations occur in either the concentration or affinity of the serum binding proteins. Therefore, in a variety of normal (pregnancy)<sup>4</sup> and abnormal (Familial Dysalbuminemic Hyperthyroxinemia, FDH)<sup>5-7</sup> states, or as a result of the administration of certain drugs (e.g. furosemide<sup>8,9</sup> and fenclofenac<sup>10-12</sup>), the target tissues are assured of receiving the required amount of hormone. Free T<sub>4</sub> values may, therefore, provide the best indication of thyroid dysfunction, since free T<sub>4</sub> is less sensitive to changes in the serum binding proteins.

Historically, the diagnosis of thyroid function has involved performing a total T<sub>4</sub> assay<sup>13,14</sup> in addition<sup>15</sup> to a Thyroxine Uptake (TU) assay of the same sample. The mathematical combination of these two assays produces a Free Thyroxine Index (FTI) which provides an indirect proportional estimate for free T<sub>4</sub>.<sup>16</sup>

Alternatively, direct assays have been developed using equilibrium dialysis,<sup>17,18</sup> ultrafiltration,<sup>19,20</sup> RIA,<sup>21</sup> and solid-phase EIA technology<sup>22</sup> to measure free T<sub>4</sub>. In these methods, separation of free and bound tracer is achieved either with a membrane, or by binding free T<sub>4</sub> to a solid phase antibody. This extraction step removes an amount of T<sub>4</sub> which is proportional to the original amount of free T<sub>4</sub> present in the patient sample. Provided that the extracted T<sub>4</sub> is less than approximately 5% of the T<sub>4</sub> in the sample, a true estimation of the free T<sub>4</sub> content can be obtained.

The ARCHITECT Free T<sub>4</sub> assay is to be used as an aid in the assessment of thyroid status.

## BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT Free T<sub>4</sub> assay is a two-step immunoassay to determine the presence of free thyroxine (Free T<sub>4</sub>) in human serum and plasma using Chemiluminescent Microparticle Immunoassay (CMIA) technology with flexible assay protocols, referred to as Chemiflex.

In the first step, sample and anti-T<sub>4</sub> coated paramagnetic microparticles are combined. Free T<sub>4</sub> (unbound) present in the sample binds to the anti-T<sub>4</sub> coated microparticles. After washing, T<sub>3</sub> acridinium labeled conjugate is added in the second step.

Pre-Trigger and Trigger Solutions are then added to the reaction mixture; the resulting chemiluminescent reaction is measured as relative light units (RLUs). An inverse relationship exists between the amount of Free T<sub>4</sub> in the sample and the RLUs detected by the ARCHITECT *i* optical system.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

## REAGENTS

### Reagent Kit, 100 Tests/500 Tests

**NOTE:** Some kit sizes are not available in all countries or for use on all ARCHITECT *i* Systems. Please contact your local distributor.

### ARCHITECT Free T<sub>4</sub> Reagent Kit (7K65)

- **MICROPARTICLES** 1 or 4 Bottle(s) (6.6 mL/27.0 mL) anti-T<sub>4</sub> (sheep) coated Microparticles in TRIS buffer with sheep IgG stabilizers. Minimum Concentration: 0.08% solids. Preservative: antimicrobial agent.
- **CONJUGATE** 1 or 4 Bottle(s) (5.9 mL/26.3 mL) T<sub>3</sub> acridinium-labeled Conjugate in MES buffer with NaCl and Triton X-100 stabilizers. Minimum concentration: 0.2 ng/mL. Preservative: ProClin.

## Other Reagents

### ARCHITECT *i* Pre-Trigger Solution

- **PRE-TRIGGER SOLUTION** Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

### ARCHITECT *i* Trigger Solution

- **TRIGGER SOLUTION** Trigger Solution containing 0.35N sodium hydroxide.

### ARCHITECT *i* Wash Buffer

**NOTE:** Bottle and volume varies based on order.

- **WASH BUFFER** Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

## WARNINGS AND PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use.
- Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

## Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens<sup>23</sup>. Biosafety Level 2<sup>24</sup> or other appropriate biosafety practices<sup>25,26</sup> should be used for materials that contain or are suspected of containing infectious agents.

The following warnings and precautions apply to this component:

- **Conjugate**



**WARNING:** H317 Contains methylisothiazolones. May cause an allergic skin reaction.

### Prevention

P261 Avoid breathing mist / vapours / spray.  
P272 Contaminated work clothing should not be allowed out of the workplace.  
P280 Wear protective gloves / protective clothing / eye protection.

### Response

P302+P352 IF ON SKIN: Wash with plenty of soap and water.  
P333+P313 If skin irritation or rash occurs: Get medical advice / attention.  
P363 Wash contaminated clothing before use.

This material and its container must be disposed of in a safe way.

- This product (the microparticles) contain sodium azide. Contact with acids liberates very toxic gas. This material and its container must be disposed of in a safe way.
- For information on the safe disposal of sodium azide and a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

## Handling Precautions

- Do not use reagent kits beyond the expiration date.
- **Do not mix reagents from different reagent kits.**
- Prior to loading the ARCHITECT Free T<sub>4</sub> Reagent Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**
- To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
- Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
- Over time, residual liquids may dry on the septum surface. These are typically dried salts which have no effect on assay efficacy.
- For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

## Storage Instructions

- 
- 2°C The ARCHITECT Free T<sub>4</sub> Reagent Kit must be stored at 2-8°C and may be used immediately after removal from 2-8°C storage.
  - When stored and handled as directed, reagents are stable until the expiration date.
  - The ARCHITECT Free T<sub>4</sub> Reagent Kit may be stored on board the ARCHITECT *i* System for a maximum of 30 days. After 30 days, the reagent kit must be discarded. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.
  - Reagents may be stored on or off the ARCHITECT *i* System. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** After reagents are removed from the system, you must initiate a scan to update the onboard stability timer.

## Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results may be invalid and may require retesting. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

- The ARCHITECT Free T<sub>4</sub> assay file must be installed on the ARCHITECT *i* System from the ARCHITECT *i* Assay CD-ROM prior to performing the assay. For detailed instructions on assay file installation and on viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.
- For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.
- For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.
- The default result unit for the ARCHITECT Free T<sub>4</sub> assay is ng/dL. An alternate result unit, pmol/L, may be selected for reporting results by editing assay parameter "Result concentration units", to pmol/L. The conversion factor used by the system is 12.87.

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

- Human serum (including serum collected in serum separator tubes) or plasma collected in sodium heparin, lithium heparin (including lithium heparin plasma separator tubes), or potassium EDTA anticoagulant tubes may be used in the ARCHITECT Free T<sub>4</sub> assay. Other anticoagulants have not been validated for use with the ARCHITECT Free T<sub>4</sub> assay. Follow the manufacturer's processing instructions for serum or plasma collection tubes.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.
- The ARCHITECT *i* System does not provide the capability to verify the specimen type. It is the responsibility of the operator to verify the correct specimen types are used in the ARCHITECT Free T<sub>4</sub> assay.
- Use caution when handling patient specimens to prevent cross contamination. Use of disposable pipettes or pipette tips is recommended.
- Do not use heat-inactivated specimens.
- For optimal results, inspect all samples for bubbles. Remove bubbles with an applicator stick prior to analysis. Use a new applicator stick for each sample to prevent cross contamination.
- For optimal results, serum and plasma specimens should be free of fibrin, red blood cells or other particulate matter.
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.

- If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator, plasma separator or red blood cells. Follow the manufacturer's processing instructions for serum or plasma collection tubes if a removal time of less than 24 hours is specified. Specimens may be stored for up to 6 days at 2-8°C prior to being tested. If testing will be delayed more than 6 days specimens should be frozen at -10°C or colder. Specimens stored frozen at -10°C or colder for 6 days showed no performance difference.
- Multiple freeze-thaw cycles of specimens should be avoided. Specimens must be mixed THOROUGHLY after thawing, by LOW speed vortexing or by gently inverting, and centrifuged prior to use to remove red blood cells or particulate matter to ensure consistency in the results.
- When shipped, specimens must be packaged and labeled in compliance with applicable state, federal and international regulations covering the transport of clinical specimens and infectious substances. Prior to shipment, it is recommended that specimens be removed from the clot, serum separator or red blood cells.

## PROCEDURE

### Materials Provided:

- 7K65 ARCHITECT Free T<sub>4</sub> Reagent Kit

### Materials Required but not Provided:

- ARCHITECT *i* System
- ARCHITECT *i* **ASSAY CD-ROM**
- 7K65-01 ARCHITECT Free T<sub>4</sub> Calibrators
- 7K65-10 ARCHITECT Free T<sub>4</sub> Controls
- ARCHITECT *i* **PRE-TRIGGER SOLUTION**
- ARCHITECT *i* **TRIGGER SOLUTION**
- ARCHITECT *i* **WASH BUFFER**
- ARCHITECT *i* **REACTION VESSELS**
- ARCHITECT *i* **SAMPLE CUPS**
- ARCHITECT *i* **SEPTUM**
- ARCHITECT *i* **REPLACEMENT CAPS**
- For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

### Assay Procedure

- Before loading the ARCHITECT Free T<sub>4</sub> Reagent Kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that have settled during shipment:
  - Invert the microparticle bottle 30 times.
  - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
  - Once the microparticles have been resuspended, remove and discard the cap. Wearing clean gloves, remove a septum from the bag. Carefully snap the septum onto the top of the bottle.
- **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
- Order tests.
- Load the ARCHITECT Free T<sub>4</sub> Reagent Kit on the ARCHITECT *i* System. Verify that all necessary assay reagents are present. Ensure that septums are present on all reagent bottles.
- The minimum sample cup volume is calculated by the system and is printed on the Orderlist report. No more than 10 replicates may be sampled from the same sample cup. To minimize the effects of evaporation verify adequate sample cup volume is present prior to running the test.
  - Priority: 95 µL for the first Free T<sub>4</sub> test plus 45 µL for each additional Free T<sub>4</sub> test from the same sample cup.
  - ≤ 3 hours on board: 150 µL for the first Free T<sub>4</sub> test plus 45 µL for each additional Free T<sub>4</sub> test from the same sample cup.
  - > 3 hours on board: additional sample volume is required. Refer to the ARCHITECT System Operations Manual, Section 5 for information on sample evaporation and volumes.

- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- ARCHITECT Free T<sub>4</sub> Calibrators and Controls should be mixed by gentle inversion prior to use.
- To obtain the recommended 150 µL volume requirements for the ARCHITECT Free T<sub>4</sub> Calibrators and Controls, hold the bottles **vertically** and dispense 4 drops of each calibrator or 4 drops of each control into each respective sample cup.
- Load samples
  - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN. The ARCHITECT *i* System performs the following function:
  - Moves the sample to the aspiration point
  - Loads a reaction vessel (RV) into the process path
  - Aspirates and transfers sample into the RV
  - Advances the RV one position and transfers microparticles into the RV
  - Mixes, incubates and washes the reaction mixture
  - Adds conjugate to the RV
  - Mixes, incubates and washes the reaction mixture
  - Adds Pre-Trigger and Trigger Solutions
  - Measures chemiluminescent emission to determine the quantity of Free T<sub>4</sub> in the sample
  - Aspirates contents of RV to liquid waste and unloads RV to solid waste
  - Calculates the result
- For information on ordering patient specimens and controls and general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- For optimal performance, it is important to follow the routine maintenance procedures defined in the ARCHITECT System Operations Manual, Section 9. If your laboratory requires more frequent maintenance, follow those procedures.

**Specimen Dilution Procedures**

Samples cannot be diluted for Free T<sub>4</sub> determinations. Samples which read > 6.00 ng/dL should be reported as such.

**Calibration**

- To perform an ARCHITECT Free T<sub>4</sub> calibration, test Calibrators 1 and 2 in duplicate. A single sample of all levels of Free T<sub>4</sub> controls must be tested to evaluate the assay calibration. Ensure that assay control values are within the concentration ranges specified in the package insert. Calibrators should be priority loaded.
- Calibrator Range: 0.0 - 6.0 ng/dL.
- Once an ARCHITECT Free T<sub>4</sub> calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
  - A reagent kit with a new lot number is used
  - Controls are out of range
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

**QUALITY CONTROL PROCEDURES**

The recommended control requirement for the ARCHITECT Free T<sub>4</sub> assay is a single sample of all control levels tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures. Ensure that assay control values are within the concentration ranges specified in the package insert.

**Verification of Assay Claims**

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B. The ARCHITECT Free T<sub>4</sub> assay belongs to method group 6.

**RESULTS**

The ARCHITECT Free T<sub>4</sub> assay utilizes a 4 Parameter Logistic Curve Fit data reduction method (4PLC, Y weighted) to generate a calibration curve.

**Alternate Result Units**

- The default result unit for the ARCHITECT Free T<sub>4</sub> assay is ng/dL. When the alternate result unit, pmol/L, is selected, the conversion factor used by the system is 12.87.
- Conversion Formula:  
(Concentration in ng/dL) x (12.87) = Concentration in pmol/L.

**Flags**

- Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

**LIMITATIONS OF THE PROCEDURE**

- For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other thyroid tests, clinical impressions, etc.
- If the Free T<sub>4</sub> results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- Performance of this test has not been established with neonatal specimens.

**EXPECTED VALUES**

A normal range of 0.70 ng/dL to 1.48 ng/dL (central 99% interval) was obtained by testing serum specimens from 411 individuals determined as normal by AxSYM Ultrasensitive hTSH II and AxSYM Free T<sub>4</sub> assays. It is recommended that each laboratory establish its own normal range, which may be unique to the population it serves depending upon geographical, patient, dietary, or environmental factors.

**SPECIFIC PERFORMANCE CHARACTERISTICS**

**Precision**

The ARCHITECT Free T<sub>4</sub> assay is designed to have a precision of ≤ 10% (total CV) for concentrations in the range of the low control (0.65 ng/dL), medium control (1.2 ng/dL) and high control (2.8 ng/dL). A study based on guidance from Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) document EP5-A<sup>27</sup> was performed for the ARCHITECT Free T<sub>4</sub> assay. A three member processed human serum based panel was assayed, using two lots of reagents, in replicates of two at two separate times per day for 20 testing days. Data from this study are summarized in the following table.\*

| Panel Member | Reagent Lot | Instrument | n  | Mean Conc. Value (ng/dL) | Within Run SD | %CV | Total SD | %CV |
|--------------|-------------|------------|----|--------------------------|---------------|-----|----------|-----|
| 1            | 1           | 1          | 80 | 0.69                     | 0.021         | 3.0 | 0.032    | 4.7 |
| 1            | 1           | 2          | 80 | 0.67                     | 0.036         | 5.3 | 0.041    | 6.1 |
| 1            | 2           | 1          | 80 | 0.70                     | 0.021         | 3.0 | 0.055    | 7.8 |
| 1            | 2           | 2          | 80 | 0.72                     | 0.027         | 3.7 | 0.043    | 6.0 |
| 2            | 1           | 1          | 80 | 1.26                     | 0.048         | 3.8 | 0.061    | 4.8 |
| 2            | 1           | 2          | 80 | 1.22                     | 0.029         | 2.3 | 0.044    | 3.6 |
| 2            | 2           | 1          | 80 | 1.25                     | 0.029         | 2.3 | 0.066    | 5.2 |
| 2            | 2           | 2          | 80 | 1.27                     | 0.033         | 2.6 | 0.048    | 3.8 |
| 3            | 1           | 1          | 80 | 2.94                     | 0.084         | 2.8 | 0.148    | 5.1 |
| 3            | 1           | 2          | 80 | 2.87                     | 0.097         | 3.4 | 0.151    | 5.3 |
| 3            | 2           | 1          | 80 | 3.03                     | 0.098         | 3.3 | 0.191    | 6.3 |
| 3            | 2           | 2          | 80 | 3.00                     | 0.088         | 2.9 | 0.134    | 4.5 |

\* Representative data; results in individual laboratories may vary from these data.

**Analytical Sensitivity**

The ARCHITECT Free T<sub>4</sub> assay is designed to have an analytical sensitivity of ≤ 0.4 ng/dL.

Analytical sensitivity is defined as the concentration calculated as the mean plus two standard deviations of replicates of the ARCHITECT Free T<sub>4</sub> MasterCheck Level 0 (0.0 ng/dL). The analytical sensitivity (low-linearity) is defined in the ARCHITECT Free T<sub>4</sub> assay parameters as 0.4 ng/dL.

### Analytical Specificity

The ARCHITECT Free T<sub>4</sub> assay is designed to have a mean analytical specificity of  $\leq 0.0035\%$  cross reactivity with triiodothyronine (T<sub>3</sub>) at a concentration of 12,000 ng/dL in a sample containing 0.5 ng/dL of Free T<sub>4</sub>.

### Interference

The ARCHITECT Free T<sub>4</sub> assay is designed to have a mean potential interference from hemoglobin, bilirubin, triglycerides, and protein of  $< 10\%$  at the levels indicated below.

- Hemoglobin -  $\leq 500$  mg/dL
- Bilirubin -  $\leq 20$  mg/dL
- Triglycerides -  $\leq 3000$  mg/dL
- Protein -  $\leq 12$  g/dL

### Accuracy by Correlation

The ARCHITECT Free T<sub>4</sub> assay is designed to have a slope of  $1.00 \pm 0.20$  and a correlation coefficient (r) of  $\geq 0.90$  when compared to the AxSYM Free T<sub>4</sub> assay.

A study was performed where specimens were tested using the ARCHITECT Free T<sub>4</sub> assay and AxSYM Free T<sub>4</sub> assay. Data from this study were analyzed using least squares and Passing Bablok<sup>28</sup> regression methods and are summarized in the following table.\*

Abbott ARCHITECT Free T<sub>4</sub> vs. Abbott AxSYM Free T<sub>4</sub>

| Method              | Number of Specimens | Intercept | Slope | Correlation Coefficient |
|---------------------|---------------------|-----------|-------|-------------------------|
| Least Squares       |                     |           |       |                         |
| Linear Regression   | 675                 | 0.03      | 0.96  | 0.953                   |
| Passing-Bablok      |                     |           |       |                         |
| Linear Regression** | 675                 | -0.02     | 1.00  | 0.953                   |

In this evaluation, serum specimens tested ranged from 0.52 ng/dL to 3.88 ng/dL with the ARCHITECT Free T<sub>4</sub> assay and from 0.46 ng/dL to 4.14 ng/dL with the AxSYM Free T<sub>4</sub> assay.

\* Representative data; variables such as differences in sampling size and sample population may impact the correlation of the assay; therefore, results in individual laboratories may vary from these data.

\*\* A linear regression method with no special assumptions regarding the distribution of the samples and the measurement errors.<sup>28</sup>

### BIBLIOGRAPHY

1. Oppenheimer JH. Role of Plasma Proteins in the Binding, Distribution and Metabolism of the Thyroid Hormones. *N Engl J Med* 1968; 278:1153-1162.
2. Woeber KA, Ingbar SH. The Contribution of Thyroxine-Binding Prealbumin to the Binding of Thyroxine in Human Serum, as Assessed by Immunoabsorption. *J Clin Invest* 1968;47:1710-1721.
3. Tabachnick M, Giorgio NA Jr. Thyroxine – Protein Interactions *Arch Biochem Biophys* 1964;105:563-569.
4. DeGroot LJ, Larsen PR, Refetoff S, Stanbury JB. Transport of Thyroid Hormone and Cell Uptake. In: *The Thyroid and Its Diseases*. New York: Wiley and Sons, 1984;62-65.
5. Stockigt JR, Topliss DJ, Barlow JW, White EL, Hurley DM, Taft P. Familial Euthyroid Thyroxine Excess: An Appropriate Response to Abnormal Thyroxine Binding Associated with Albumin. *J Clin Endocrinol Metab* 1981;53:353-359.
6. Ruiz M, Rajatanavin R, Young RA, Taylor C, Brown R, Braverman LE, et al. Familial Dysalbuminemic Hyperthyroxinemia. *N Engl J Med* 1982;306:635-639.
7. Hennemann G, Docter R, Krenning EP, Bos G, Otten M, Visser TJ. Raised Total Thyroxine and Free Thyroxine Index but Normal Free Thyroxine. *Lancet* 1979;1:639-642.
8. Stockigt JR, Lim CF, Barlow JW, Stevens V, Topliss DJ, Wynne KN. High Concentrations of Furosemide Inhibit Serum Binding of Thyroxine. *J Clin Endocrinol Metab* 1984;59:62-66.
9. Stockigt JR, Lim CF, Barlow JW, Stevens V, Wynne KN, Mohr VS, Topliss DJ, et al. Interaction of Furosemide with Serum Thyroxine-Binding Sites: *In Vivo* and *In Vitro* Studies and Comparison with Other Inhibitors. *J Clin Endocrinol Metab* 1985;60:1025-1031.
10. Ratcliffe WA, Hazelton RA, Thomson JA, Ratcliffe JG. The Effect of Fenclufenac On Thyroid Function Tests *In Vivo* and *In Vitro*. *Clin Endocrinol* 1980;13:569-575.

11. Kurtz AB, Capper SJ, Clifford J, Humphrey MJ, Lukinac L. The Effect of Fenclufenac On Thyroid Function. *Clin Endocrinol* 1981;15:117-124.
12. Pearson DWM, Ratcliffe WA, Thomson JA, Ratcliffe JG. Biochemical and Clinical Effects of Fenclufenac In Thyrotoxicosis. *Clin Endocrinol* 1982;16:369-373.
13. Larsen PR, Dockalova J, Sipula D, Wu FM. Immunoassay of Thyroxine in Unextracted Human Serum. *J Clin Endocrinol Metab* 1973; 37:177-182.
14. Chopra IJ. A Radioimmunoassay for Measurement of Thyroxine in Unextracted Human Serum. *J Clin Endocrinol Metab* 1972; 34:938-947.
15. Refetoff S, Hayen SR, Selenkow HA. Estimation of the T<sub>4</sub> Binding Capacity of Serum TBG and TBPA by a Single T<sub>4</sub> Load Ion Exchange Resin Method. *J Nucl Med* 1972;13:2-12.
16. Clark F, Horn DB. Assessment of Thyroid function by the Combined Use of the Serum Protein-Bound Iodine and Resin Uptake of 131 I-Triiodothyronine. *J Clin Endocrinol Metab* 1965;25:39-45.
17. Ellis SM, Ekins RP. The Radioimmunoassay of Serum Free Triiodothyronine and Thyroxine. In: Pasternak CA, editor. *Radioimmunoassay in Clinical Biochemistry*. London: Heyden, 1975;25:187-194.
18. Nelson JC, Tomei RT. Direct Determination of Free Thyroxine in Undiluted Serum by Equilibrium Dialysis/Radioimmunoassay in Serum. *Clin Chem* 1988;34:1737-1744.
19. Sophianopoulos J, Jerkunica I, Lee CN, Sgoutas D. An Improved Ultrafiltration Method for Free Thyroxine and Triiodothyronine in Serum. *Clin Chem* 1980;26:159-162.
20. Shannon N, Woolf PD. Determination of Free Thyroxine in serum by Ultrafiltration: Validation of a Method and Preliminary Results. *Clin Chem* 1984;30:1770-1773.
21. Bayer MF, McDougall IR. Radioimmunoassay of Free Thyroxine in Serum: Comparison with Clinical Findings and Results of Conventional Thyroid-Function Tests. *Clin Chem* 1980;26:1186-1192.
22. Sturgess ML, Weeks I, Evans PJ, Mpoko CN, Laing I, Woodhead JS. An Immunochemiluminometric Assay for Serum Free Thyroxine. *Clin Endocrinol* 1987;27:383-393.
23. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
24. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; January 2007.
25. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
26. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline - Third Edition*. CLSI Document M29-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
27. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices – Approved Guideline. NCCLS Document EP5-A. Wayne, PA: NCCLS, 1999.
28. Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. *J Clin Chem. Clin Biochem*. 1983;21:709-720.

ARCHITECT, AxSYM, Chemiflex and MasterCheck are trademarks of Abbott Laboratories in various jurisdictions.

All trademarks are property of their respective owners.



Abbott Ireland  
Diagnostics Division  
Lisnamuck, Longford  
Co. Longford  
Ireland  
+353-43-3331000



Distributed by Abbott Laboratories  
Abbott Park, IL 60064 USA  
and  
ABBOTT 65205 Wiesbaden, Germany

January 2010  
© 2006, 2010 Abbott Laboratories