

ARCHITECT / AEROSET

LIPASE

This package insert contains information to run the Lipase assay on the ARCHITECT c Systems and the AEROSET System.






NOTE: Changes Highlighted

NOTE: This package insert must be read carefully prior to product use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Customer Support

United States: 1-877-4ABBOTT
 Canada: 1-800-387-8378 (English speaking customers)
 1-800-465-2675 (French speaking customers)
 International: Call your local Abbott representative

Symbols in Product Labeling

CAL	Calibrator	REF	Catalog number/List number
CONC	Concentration	SN	Serial number
EC REP	Authorized Representative in the European Community		Consult instructions for use
INGRED	Ingredients		Manufacturer
IVD	In vitro diagnostic medical device		Temperature limitation
LOT	Batch code/Lot number		Use by/Expiration date
R1	Reagent 1		CAUTION. Consult accompanying documents.
R1A	Reagent 1A		
R2	Reagent 2		

NAME

LIPASE

INTENDED USE

The Lipase assay is used for the quantitation of lipase in human serum or plasma.

SUMMARY AND EXPLANATION OF TEST

Pancreatic lipase in serum and plasma is closely associated with pancreatic diseases. The activity of this enzyme has been measured as an important marker for diagnosing pancreatic diseases and the associated monitoring of therapeutic effects. Pancreatic lipase test kits currently available include a turbidimetric method using triglyceride as substrate and a colorimetric method using synthetic substrates. These methods, however: 1) lack precision near the normal level; 2) exhibit poor reproducibility; and 3) are affected by other enzymes such as esterases.

The enzymatic color rate assay uses a clear substrate solution of 1,2-diglyceride, which is a 'natural' substrate. The assay is highly sensitive and specific for pancreatic lipase, using colipase and deoxycholate as activators.

PRINCIPLES OF PROCEDURE

Lipase acts on a natural substrate, 1,2-diglyceride, to liberate 2-monoglyceride. This is hydrolyzed by monoglyceride lipase into glycerol and free fatty acid. Glycerol kinase acts on glycerol to form glycerol-3-phosphate which in turn acted on by glycerol-3-phosphate oxidase to generate hydrogen peroxide. Peroxidase converts the hydrogen peroxide, 4-aminoantipyrine, and *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-*m*-toluidine (TOOS) into a quinone dye. The rate of formation of the dye, measured as an increase in absorbance at 548 nm, is proportional to the lipase concentration in the sample.

Methodology: Quinone Dye

REAGENTS

Reagent Kit

[REF] 7D80 Lipase is supplied as a two-reagent kit which contains:

[R1] 5 x 29 mL

[R1A] 5 x 30 mL
Contains human serum albumin.

[R2] 5 x 14 mL

Estimated tests per kit: 778

Calculation is based on the minimum reagent fill volume per kit.

Reactive Ingredients	Concentration
[R1] Cholic Acid	5.34 mmol/L
[R1A] 1,2-Diglyceride	1.1 mmol/L
Monoglyceride Lipase	≥ 0.86 U/mL
Glycerol Kinase	≥ 1.34 U/mL
Glycerol-3-Phosphate Oxidase	≥ 40.0 U/mL
Peroxidase	≥ 1.34 U/mL
Colipase	≥ 40.0 U/mL
TOOS	0.068%
ATP	0.66 mmol/L
[R2] Deoxycholate	36.0 mmol/L
4-Aminoantipyrine	0.12%

REAGENT HANDLING AND STORAGE

Reagent Handling

Remove air bubbles, if present in the reagent cartridge, with a new applicator stick. Alternatively, allow the reagent to sit at the appropriate storage temperature to allow the bubbles to dissipate. To minimize volume depletion, do not use a transfer pipette to remove the bubbles.

CAUTION: Reagent bubbles may interfere with proper detection of reagent level in the cartridge, causing insufficient reagent aspiration which could impact results.

Reagent Storage


Unopened reagents are stable until the expiration date when stored at 2 to 8°C.

Reagent stability is 11 days if the reagent is uncapped and onboard.

REAGENT HANDLING AND STORAGE (Continued)

Instructions for Use

1. Prepare the Working Reagent by adding the contents of the large Lipase [R1] to the smaller Lipase [R1A].



2. Replace the [R1A] reagent stopper and mix by gentle inversion to achieve complete dissolution.
3. Return the prepared Working Reagent to the [R1] cartridge and mix again by gentle inversion.
Remove air bubbles, if present in the reagent cartridge, with a new applicator stick.
4. Place the [R1] cartridge in Reagent Supply Center 1.
5. Place the [R2] cartridge in Reagent Supply Center 2.

WARNINGS AND PRECAUTIONS

Precautions for Users

1. For in vitro diagnostic use.
2. Do not use components beyond the expiration date.
3. Do not mix materials from different kit lot numbers.
4. Do not mix reagents prepared at different times.
5. Do not reuse the reagent containers, bottles, caps, or plugs due to the risk of contamination and the potential to compromise reagent performance.
6. **CAUTION:** This product contains human sourced and/or potentially infectious components. Refer to the REAGENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.
The human serum albumin used in Reagent 1A has been tested and found to be nonreactive for HBsAg, anti-HCV, and anti-HIV-1/HIV-2.

7. [R1] and [R2] contain sodium azide. Contact with acids liberates very toxic gas. This material and its container must be disposed of in a safe way.

For reagents not classified as dangerous per European Directive 1999/45/EC as amended, safety data sheet available for professional user on request.

SPECIMEN COLLECTION AND HANDLING

Suitable Specimens

Serum and plasma are acceptable specimens.

- **Serum:** Use serum collected by standard venipuncture techniques into glass or plastic tubes with or without gel barriers. Ensure complete clot formation has taken place prior to centrifugation. When processing samples, separate serum from blood cells or gel according to the specimen collection tube manufacturer's instructions.
Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may take longer to complete their clotting processes. Fibrin clots may subsequently form in these sera and the clots could cause erroneous test results.
- **Plasma:** Use plasma collected by standard venipuncture techniques into glass or plastic tubes. Acceptable anticoagulants are lithium heparin (with or without gel barrier) and sodium heparin. Ensure centrifugation is adequate to remove platelets. When processing samples, separate plasma from blood cells or gel according to the specimen collection tube manufacturer's instructions.

For total sample volume requirements, refer to the instrument-specific ASSAY PARAMETERS section of this package insert and Section 5 of the instrument-specific operations manual.

SPECIMEN COLLECTION AND HANDLING (Continued)

Specimen Storage

Serum/Plasma

Temperature	Maximum Storage	Bibliographic Reference
20 to 25°C	7 days	5
2 to 8°C	7 days	5, 6
-20°C	1 year	5

Guder et al.⁵ suggest storage of frozen specimens at -20°C for no longer than the time interval cited above. However, limitations of laboratory equipment make it necessary in practice for clinical laboratories to establish a range around -20°C for specimen storage. This temperature range may be established from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

NOTE: Stored specimens must be inspected for particulates. If present, mix and centrifuge the specimen to remove particulates prior to testing.

PROCEDURE

Materials Provided

[REF] 7D80 Lipase Reagent Kit

Materials Required but not Provided

- [REF] 3E16 Lipase Calibrator, [CAL] 2 x 3 mL
- Control Material
- Saline (0.85% to 0.90% NaCl) for specimens that require dilution
- [REF] 2J94 Detergent B for ARCHITECT c Systems only

Assay Procedure

For a detailed description of how to run an assay, refer to *Section 5* of the instrument-specific operations manual.

Specimen Dilution Procedures

The ARCHITECT c Systems and the AEROSET System have automatic dilution features; refer to *Section 2* of the instrument-specific operations manual for additional information.

Serum and plasma: Specimens with lipase values exceeding 1,200 U/L are flagged and may be diluted using the Automated Dilution Protocol or the Manual Dilution Procedure.

Automated Dilution Protocol

If using the Automated Dilution Protocol, the system performs a dilution of the specimen and automatically corrects the enzyme activity value by multiplying the result by the appropriate dilution factor. To set up the automatic dilution feature, refer to *Section 2* of the instrument-specific operations manual for additional information.

Manual Dilution Procedure

Manual dilutions should be performed as follows:

- Use saline (0.85% to 0.90% NaCl) to dilute the sample.
- The operator must enter the dilution factor in the patient or control order screen. The system uses this dilution factor to automatically correct the enzyme activity value by multiplying the result by the entered factor.
- If the operator does not enter the dilution factor, the result must be multiplied by the appropriate dilution factor before reporting the result.

NOTE: If a diluted sample result is flagged indicating it is less than the linear low limit, do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to *Section 5* of the instrument-specific operations manual.

CALIBRATION

Calibration is stable for approximately 11 days (264 hours) and is required with each change in reagent lot number. Verify calibration with at least two levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

For a detailed description of how to calibrate an assay, refer to *Section 6* of the instrument-specific operations manual.

For information on calibrator standardization, refer to the Lipase Calibrator package insert.

QUALITY CONTROL

The following is the recommendation of Abbott Laboratories for quality control. As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls (normal and abnormal) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, patient values may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

RESULTS

Refer to the instrument-specific operations manual for information on results calculations.

- **ARCHITECT System Operations Manual—Appendix C**
- **AEROSET System Operations Manual—Appendix A**

Representative performance data are given in the EXPECTED VALUES and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert. Results obtained in individual laboratories may vary.

LIMITATIONS OF THE PROCEDURE

Refer to the SPECIMEN COLLECTION AND HANDLING and PERFORMANCE CHARACTERISTICS sections of this package insert.

For the ARCHITECT c Systems only: Detergent B solution must be onboard the system.

For the AEROSET System only: Lipase must be run on a separate line from Calcium, Cholesterol, Carbon Dioxide, MULTIGENT Direct LDL, Triglyceride, and Uric Acid. Refer to *Section 2* of the **AEROSET System Operations Manual** for additional information.

EXPECTED VALUES

Reference Range

Serum⁷

	Range (U/L)
Adult	8 to 78

A study was conducted using 133 serum samples from volunteers. Data were analyzed as described by Clinical and Laboratory Standards Institute (CLSI) protocol NCCLS C28-A.⁸ From this study, 95% of all specimens fell within 8 to 78 U/L, with samples ranging from 7 to 90 U/L.

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

SPECIFIC PERFORMANCE CHARACTERISTICS

Linearity

Lipase is linear up to 1,200 U/L. Linearity was verified using CLSI protocol NCCLS EP6-P.⁹

Limit of Detection (LOD)

The LOD for Lipase is 1.6 U/L. The LOD is the mean concentration of an analyte-free sample + 2 SD, where SD = the pooled, within-run standard deviation of the analyte-free sample.

Limit of Quantitation (LOQ)

The LOQ for Lipase is 3.1 U/L. The LOQ is the analyte concentration at which the CV = 20%.

SPECIFIC PERFORMANCE CHARACTERISTICS

(Continued)

Interfering Substances

Interference studies were conducted using CLSI protocol NCCLS EP7-P.¹⁰ Interference effects were assessed by Dose Response and Paired Difference methods, at the medical decision level of the analyte.

Interfering Substance	Interferent Concentration	N	Target (U/L)	Observed (% of Target)
Bilirubin	15 mg/dL (257 µmol/L)	4	212.3	94.7
	30 mg/dL (513 µmol/L)	4	212.3	86.4
Hemoglobin	1,000 mg/dL (10.0 g/L)	4	215.7	103.1
	2,000 mg/dL (20.0 g/L)	4	215.7	99.1
Intralipid	750 mg/dL (7.5 g/L)	4	198.1	94.1
	1,000 mg/dL (10.0 g/L)	4	198.1	93.9

Bilirubin solutions at the above concentrations were prepared by addition of a bilirubin stock to human serum pools. Hemoglobin solutions at the above concentrations were prepared by addition of hemolysate to human serum pools. Intralipid solutions at the above concentrations were prepared by addition of Intralipid to human serum pools.

Interferences from medications or endogenous substances may affect results.¹¹

Precision

The imprecision of the Lipase assay is $\leq 7.5\%$ Total CV. Representative data from studies using CLSI protocol NCCLS EP5-T2¹² are summarized below.

Control		Level 1	Level 2
N		80	80
Mean (U/L)		42.7	78.6
Within Run	SD	0.31	1.43
	%CV	0.7	1.8
Between Run	SD	1.13	1.88
	%CV	2.6	2.4
Between Day	SD	0.91	2.62
	%CV	2.1	3.3
Total	SD	1.48	3.53
	%CV	3.5	4.5

Method Comparison

Correlation studies were performed using CLSI protocol NCCLS EP9-A.¹³

Serum results from the Lipase assay on the AEROSSET System were compared with those from a commercially available enzymatic methodology.

Serum results from the Lipase assay on an ARCHITECT cSystem were compared with the Lipase assay on the AEROSSET System.

	AEROSSET vs. Comparative Method	ARCHITECT vs. AEROSSET
N	74	80
Y - Intercept	-1.458	2.891
Correlation Coefficient	1.000	0.997
Slope	0.936	1.032
Range (U/L)*	10.3 to 619.1	3.10 to 1,146.80

*AEROSSET

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030, *Bloodborne Pathogens*.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, January 2007.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Sewell DL, Bove KE, Callihan DR, et al. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition (M29-A3)*. Wayne, PA: Clinical and Laboratory Standards Institute, 2005.
5. Guder WG, da Fonseca-Wollheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2001:36–7.
6. US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
7. Data on file at Abbott Laboratories.
8. Sasse EA, Aziz KJ, Harris EK, et al. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline (C28-A)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1995.
9. Passey RB, Bee DE, Caffo A, et al. *Evaluation of the Linearity of Quantitative Analytical Methods; Proposed Guideline (EP6-P)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1986.
10. Powers DM, Boyd JC, Glick MR, et al. *Interference Testing in Clinical Chemistry; Proposed Guideline (EP7-P)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1986.
11. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3-398–3-400.
12. Kennedy JW, Carey RN, Coolen RB, et al. *Evaluation of Precision Performance of Clinical Chemistry Devices—Second Edition; Tentative Guideline (EP5-T2)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1992.
13. Kennedy JW, Carey RN, Coolen RB, et al. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline (EP9-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 1995.

TRADEMARKS

The ARCHITECT cSystem family of instruments consists of c4000, c8000, and c16000 instruments.

AEROSSET, ARCHITECT, c4000, c8000, c16000, cSystem, MULTIGENT, and SmartWash are trademarks of Abbott Laboratories.

All other trademarks, brands, product names, and trade names are the property of their respective companies.

ARCHITECT c SYSTEMS ASSAY PARAMETERS

ARCHITECT

Lipase Serum/Plasma—Conventional and SI Units

Configure assay parameters — General				
<input checked="" type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation				
Assay: Lip		Type: Photometric		Version: †
Number: 1029				
<input checked="" type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input type="radio"/> Validity checks				
Reaction mode: Rate up				
Primary		Secondary		Read times
Wavelength: 548 / 660				
Last required read: 33		Main: 22 – 33		
Absorbance range: ___-___		Flex: ___-___		
Sample blank type: None		Color Correction: ___-___		

Configure assay parameters — Reagent / Sample				
<input type="radio"/> Reaction definition <input checked="" type="radio"/> Reagent / Sample <input type="radio"/> Validity checks				
Reagent: LIP00		Reagent volume: 156		R1 R2
Diluent: Saline		Water volume: ___		
Diluent dispense mode: Type 0				
Dispense mode: Type 0		Type 0		
Dilution name	Sample	Diluted sample	Diluent	Water
STANDARD	4.0	___	___	___
___	___	___	___	___
___	___	___	___	___
Dilution factor = 1:1.00				Default dilution
				<input checked="" type="radio"/>
				<input type="radio"/>
				<input type="radio"/>

Configure assay parameters — Validity checks				
<input type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input checked="" type="radio"/> Validity checks				
Reaction check: End Subtraction				
		A		B
		Read Time: 33 – 33		31 – 31
		Calculation limits: 0.0001 – 9.9999		
Rate linearity %: ___				

Configure assay parameters — Calibration				
<input type="radio"/> General <input checked="" type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation				
Assay: Lip		Calibration method: Linear		
<input checked="" type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks				
Calibrator set: Lipase		Calibrator level: Water		Concentration: 0
Replicates: 3 [Range 1 – 3]		Cal 1: Lipase1		‡

Configure assay parameters — Volumes				
<input type="radio"/> Calibrators <input checked="" type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks				
Calibrator: Lipase				
	Calibrator level	Sample	Diluted sample	Diluent
	Blank: Water	4.0	___	___
	Cal 1: Lipase1	4.0	___	___

Configure assay parameters — Intervals				
<input type="radio"/> Calibrators <input type="radio"/> Volumes <input checked="" type="radio"/> Intervals <input type="radio"/> Validity checks				
Calibration intervals:				
Full interval: 264		(hours)		
Calibration type:				
Adjust type: None				

Configure assay parameters — Validity checks				
<input type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input checked="" type="radio"/> Validity checks				
Blank absorbance range: ___ - ___				
Span: Blank		- Blank		
Span absorbance range: ___ - ___				
Expected cal factor: 0.00				
Expected cal factor tolerance %: 0				

Configure assay parameters — SmartWash				
<input type="radio"/> General <input type="radio"/> Calibration <input checked="" type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation				
Assay: Lip				
COMPONENT	REAGENT / ASSAY	WASH	Volume	Replicates
R1	ACETM	0.5% AcidWash	345	1†
R1	LIP00	10% Detergent B	345	2
R1	All	10% Detergent B	345	2
R2	LIP00	10% Detergent B	345	2
R2	All	10% Detergent B	345	2
Cuvette	CO2	10% Detergent B	345	
Cuvette	Uric	10% Detergent B	345	
Cuvette	Trig	10% Detergent B	345	
Cuvette	Chol	10% Detergent B	345	
Cuvette	Uric-U	10% Detergent B	345	
Sample Probe		0.5% AcidWash		
*Reagent Probe SmartWashes must be configured in order listed.				

Configure assay parameters — Results				
<input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input checked="" type="radio"/> Results <input type="radio"/> Interpretation				
Assay: Lip		Assay number: 1029		
Dilution default range:		Result units: U/L		
		Low-Linearity: 4††		
		High-Linearity: 1200		
Gender and age specific ranges:				
GENDER	AGE (UNITS)	NORMAL	EXTREME	
Either	0 – 130 (Y)	8 – 78		

Configure result units	
Assay: Lip	Version: †
Result units: U/L	
Decimal places: 0	[Range 0 – 4]
Correlation factor: 1.0000	
Intercept: 0.0000	

† Due to differences in instrument systems and unit configurations, version numbers may vary.
 †† The linear low value (Low-Linearity) is LOQ rounded up to the number of decimal places defined in the decimal places parameter field.
 ‡ Refer to concentration specified on calibrator labeling or value sheet. In ARCHITECT software version 5.00 and above, these values are defined on the Configure calibrator set screen.

AEROSET SYSTEM ASSAY PARAMETERS

AEROSET

Lipase Serum/Plasma—Conventional and SI Units

Assay Configuration: Outline Page						
Assay Name	Assay #					Line
Lip	29					A-Line
Quantitative Ranges						
Min Text	Min	Panic-L	L-Reference-H	Panic-H	Max	Max Text
*	0.0*	0.0	8	78	0.0	0.0*
		4**	L-Linear Range-H		1200	
Reference Ranges*						
	Age		Male		Female	
	0 Year		0.0 - 0.0		0.0 - 0.0	
	0 Year		0.0 - 0.0		0.0 - 0.0	
	0 Year		0.0 - 0.0		0.0 - 0.0	
Qualitative Ranges						
N/A						

Assay Configuration: Base Page			
Reaction Mode	Wavelength-Prim/Sec	Read time-Main/Flex	Linearity%
RATE UP	548 / 660	22 - 33 / 0 - 0	0
Sample Blank Test	Blank Read Time	Abs Window	Abs Limits
____ (____)	0 - 0	0 - 0	0.0 - 0.0
Standard	S.Vol	DS.Vol	D.Vol
4.0	0.0	0	0
Dil 1	4.0	0.0	0
Dil 2	4.0	0.0	0
Reagent 1	Rgt Name/Pos	R.Vol	W.Vol
LIP0051 - ____*		156	0
Reagent 2	Rgt Name/Pos	R.Vol	W.Vol
LIP0052 - ____*		52	0
Reaction Check	Read Time - A/B	Range	Minimum
END SUB	33 - 33 / 31 - 31	0.0001 - 9.9999	0.0
Factor/Intercept	Decimal Places	Units	
1.0 / 0.0	0	U/L	

Assay Configuration: Calibration Page				
Calib Mode	Interval (H)			
Linear	264			
Blank/Calib Replicates	Extrapolation%	Span	Span Abs Range	
3 / 3	0	BLK - 1	0.0 - 0.0	
Sample	S.Vol	DS.Vol	D.Vol	W.Vol
BLK Water	4.0	0.0	0	0
C1 Lipase	4.0	0.0	0	0
C2	2.0	0.0	0	0
				BLK Abs Range
				0.0 - 0.0
				Cal Deviation
				0.0
				FAC Limit (%)
				10

Assay Configuration: SmartWash Page			
Rgt Probe	Reagent***	Wash	Vol
	LIP0051	AlkW	300
	LIP0052	AlkW	300
	CK00051†	Water	345
	CK00042††	Water	345
	UHDL061	AlkW	345
	ACETM41	AcidW	345
	ACETM42	AcidW	345
	DIG0051	AlkW	345
	DIG0012	AlkW	345
	AMIK941	AlkW	345
	AMIK942	AlkW	345
	VANCO51	AlkW	345
	VANCO52	AlkW	345
	TOBRA41	AlkW	345
	TOBRA42	AlkW	345
	DGTOB11	AlkW	345
	DGTOB12	AlkW	345
Cuvette	Assay Name	Wash	Vol
	—	—	—
Sample Probe	Wash		
	—		

Refer to **Assay Configuration** in *Section 2* of the **AEROSET System Operations Manual** for information regarding assay parameters.

* User defined or instrument defined.

** The linear low value (L-Linear Range) is LOQ rounded up to the number of decimal places defined in the decimal places parameter field.

*** Rgt Probe SmartWashes must be implemented if the Lipase assay is configured on the same line as the donor assay. AEROSET software allows configuration of only 16 Reagent Probe washes per assay.

† Reagent SmartWash listed is for CK **7D63-21**. For CK **7D63-20** or **7D63-30**, change this Reagent SmartWash to **CK00061**.

†† Reagent SmartWash listed is for CK **7D63-21**. For CK **7D63-20** or **7D63-30**, change this Reagent SmartWash to **CK00052**.

